roeko

GuttaFlow® bioseal

Step by step

DE	Gebrauchsinformation	4	PL	Instrukcja użycia	30
EN	Instructions for use	6	SL	Navodila za uporabo	32
FR	Mode d'emploi	8	LT	Naudojimo instrukcijos	34
IT	Istruzioni per l'uso	10	ET	Kasutusjuhend	36
ES	Instrucciones de uso	12	CS	Návod k použití	38
PT	Instruções de utilização	14	HU	Használati utasítás	40
NL	Gebruiksaanwijzing	16	SK	Instrukcie za upotreba	42
DA	Brugsanvisning	18	BG	Návod na použitie	44
SV	Bruksanvisning	20	LV	Lietošanas norādījumi	46
NO	Bruksanvisning	22	TR	Kullanma talimatı	48
FI	Käyttöohjeet	24	RU	Инструкция по применению	50
EL	Οδηγίες χρήσης	26	HR	Upute za upotrebu	52
RO	Instrucțiuni de utilizare	28	UK	Інструкція для застосування	54

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

Removal technique for root post preparation

STEP 1:

Root canal filled with GuttaFlow® bioseal and Mastercone.


STEP 2:

Procedure:

After the GuttaFlow® bioseal has set, remove the upper part of the filling using a steel or carbide round bur without exerting pressure or force. If there is sufficient gutta-percha, you can remove it with a hot instrument.

Material:

The size of the round bur must be suitable for the clinical situation and able to keep wall contact. Work can be performed with or without water cooling. Recommended speed: turbine up to 30,000 rpm, or with red contra-angle up to 100,000 rpm (work with a soft pecking motion).

STEP 3:


Procedure:

Remove the root filling peripherally to a depth of approx. 4-5 mm using a round end tapered diamond bur head, without exerting pressure or force, in such a way that the GuttaFlow® bioseal does not wrap itself round the bur.

Remove the gutta-percha carefully along the root canal in order to prevent perforation and drilling into the canal walls (via falsa). With premolars and molars adapt the drilling depth as required.

Material:

The grain size of the round end tapered diamond bur head must be coarse. Work should be performed with water cooling. Recommended speed: turbine up to approx. 40-50,000 rpm (work with a soft pecking motion).

STEP 4:

Select the post (we recommend Para-Post®) to suit the clinical situation. Prepare the root canal in accordance with the post Instructions for Use. Recommended speed: use the turbine at up to approx. 750-1,000 rpm.

The post drill should rotate clockwise continuously until it is completely removed from the tooth. This reduces the risk of the drill binding in the canal.

STEP 5:

Introduce the post to the canal and shorten the post if necessary as explained in accordance with the post Instructions for Use. Cement the post in place in accordance with the manufacturer's instructions.

Immediate post preparation for temporary posts

EN


STEP 1:

Procedure:

If there is sufficient gutta-percha you can remove the upper part with a hot instrument.

STEP 2:

Remove the unset root filling to a depth of 3-4 mm using a round end tapered diamond. Keep wall contact. No water cooling is necessary. The heat of the diamond cures the material immediately so that it can easily be removed from the canal walls.

STEP 3:

Select the post drill which most closely matches the diameter of the tapered diamond bur used in step 2. Prepare post space in accordance to the manufacturer's instructions. Insert the corresponding temporary post into the prepared post space. If necessary, shorten from the apical end. To preserve the integrity of the post space, apply temporary cement only on the margins of the crown.

Do not place any cement into the post space or on the post.

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Coltene/Whaledent GmbH + Co. KG