
3M | Elipar™ DeepCure-L

ESPE

- ① **LED Curing Light**
- ② **LED-Lichtgerät**
- ③ **Lampe à photopolymériser LED**
- ④ **Lampada fotopolimerizzatrice LED**
- ⑤ **Lámpara de fotopolimerización dental de tecnología LED**
- ⑥ **Fotoemissor LED**
- ⑦ **LED Uithardingslamp**
- ⑧ **Συσκευή φωτοπολυμερισμού τύπου LED**
- ⑨ **LED-lampa**
- ⑩ **LED-valokovetin**
- ⑪ **LED polymeriseringslampe**
- ⑫ **LED-herdelampe**

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Elipar™ DeepCure-L

LED Curing Light

Table of Contents	Page
Safety Instructions	1
Product Description	3
Fields of Application	3
Technical Data	3
Charger	3
Handpiece	4
Charger and Handpiece	4
Transport and Storage Conditions	4
Installation of the Unit	4
Factory Settings	4
Initial Steps	4
Charger	4
Light Guide/Handpiece	5
Battery Charging	5
Battery Power Level Display on Handpiece	6
Operation	6
Selection of Exposure Time	6
Activating and Deactivating the Light	6
Inserting and Removing the Light Guide from/into the Handpiece	6
Positioning the Light Guide	7
Testing of Light Intensity	7
Recommended Cure Techniques with 3M Adhesives	7
Sleep Mode	8
Audible Signals – Handpiece	8
Troubleshooting	8
Maintenance and Care	9
Care of the Handpiece	9
Reprocessing the Light Guide	9
Clean Handpiece and Glare Shield	12
Storage of the Handpiece during Extended Periods of Non-Use	12
Return of Old Electric and Electronic Equipment for Disposal	12
Customer Information	12
Warranty	12
Limitation of Liability	12
Symbol Glossary	13

Safety Instructions

PLEASE NOTE!

Prior to installation and start-up of the device, please read these instructions carefully.

As with all technical devices, the proper function and safe operation of this device depend on the user's compliance with the standard safety procedures as well as the specific safety recommendations presented in these Operating Instructions.


1. Use of the device is restricted to trained personnel in accordance with the instructions below. The manufacturer assumes no liability for any damage arising from any other or improper use of this device.
2. The charger must be accessible at all times. Do not use the charger for any use other than the charging of the Elipar™ DeepCure-L handpiece. Disconnect the handpiece from the mains by unplugging the charger from the electrical outlet.

Treating patients using the handpiece while it is still connected to the charger is prohibited for safety reasons. Light-curing is possible only if the charger has been disconnected.

3. Use only the charger "CPS008050100" or "GTM 96060-0606-1.0" (AC adapter plug) which is provided with the device. The use of any other charger can result in damage to the battery.
4. CAUTION! Avoid using this device in close proximity to other equipment or with other equipment in stacked form, as this may result in improper function. If use in the manner described before is necessary, this device and the other equipment should be observed carefully in order to ensure that they work properly.
5. CAUTION! Do not stare at source. May be harmful to the eyes. Restrict exposure to the area of the oral cavity in which clinical treatment is intended. Protect patient and user from reflection and intensive scattered light by taking the appropriate precautions, e.g., glare shields, goggles, or coverings.
6. CAUTION! As is the case for all high-intensity light-curing devices, the high light intensity is accompanied by heat generation on the exposed surface. This heat can result in irreversible damage if there is longer exposure in the proximity of the pulp or soft tissue. The exposure times given in the manufacturer's instructions must be observed exactly to avoid any such damage.

Uninterrupted exposure times of the same tooth surface in excess of 20 seconds and direct contact with oral mucosa or skin must be strictly avoided. Scientists working in this field are in agreement that the irritation caused by heat generated during light curing can be minimized by taking two simple precautions:

- Polymerization with external cooling from an air flow
 - Polymerization at intermittent intervals (e.g., 2 exposures lasting 10 seconds each instead of 1 exposure lasting 20 seconds).
7. Elipar DeepCure-L may be operated only with the supplied light guide or original 3M Elipar DeepCure-L replacement and accessory light guide.
The light guide has to be seen as an applied part. The use of other light guides may result in a reduction or increase in the light intensity. The product's warranty does not cover any damage resulting from the use of third-party light guides.
 8. Condensation resulting from the device being transferred from a cold to a warm environment may be a potential risk. Never begin operating the device until it has reached the ambient temperature.
 9. In order to avoid electric shock, do not introduce any objects into the device with the exception of replacement parts handled in accordance with the Operating Instructions.
 10. Use only genuine 3M parts when replacing defective components as directed in these Operating Instructions. The product's warranty does not cover any damage resulting from the use of third-party replacement parts.
 11. Should you have any reason to suspect the safety of the device to be compromised, the device must be taken out of operation and labeled accordingly to prevent third parties from inadvertently using a possibly defective device. Safety may be compromised, e.g., if the device malfunctions or is noticeably damaged.
 12. Keep solvents, flammable liquids, and sources of intense heat away from the device as they may damage the plastic housing of the device, the seals, or the operating buttons.
 13. Do not operate the device in the proximity of flammable mixtures.
 14. Do not allow any cleaning agents to enter the device during cleaning as they could cause an electrical short or a dangerous malfunction.
 15. Only service centers authorized by 3M Deutschland GmbH may open the device housing and repair the device.
 16. **WARNING:** this device must not be modified without the manufacturer's permission.
 17. Elipar DeepCure-L must not be used in patients, or by users, with heart pacemaker implants who have been advised to be cautious with regard to their exposure to small electrical devices.
 18. Do not use Elipar DeepCure-L in patients with a history of photobiological reactions (including individuals with urticaria solaris or erythropoietic protoporphyria) or who are currently on photosensitizing medication (including 8-methoxypsoralen or dimethylchlorotetracycline).
 19. Individuals with a history of cataract surgery may be particularly sensitive to the exposure to light and should be discouraged from Elipar DeepCure-L treatment unless adequate safety measures, such as the use of protective goggles to remove blue light, are undertaken.
 20. Individuals with a history of retinal disease should seek advice from their ophthalmologist prior to operating the device. In operating the Elipar DeepCure-L device, this group of individuals must take extreme care and comply with any and all safety precautions (including the use of suitable light-filtering safety goggles).
 21. This device has been developed and tested in accordance with the relevant EMC regulations and standards. It is in conformity with legal requirements. Since various factors such as power supply, wiring, and the ambient conditions at the place of operation can affect the EMC properties of the device, the possibility that, under unfavorable conditions, there will be EMC disruptions cannot be completely excluded. If you should notice problems in the operation of this or other devices, move the device to a different location. The EMC manufacturer's declaration and the recommended separation distances between portable and mobile RF communications equipment and the Elipar DeepCure-L unit are listed in the appendix.
 22. **CAUTION!** Portable high-frequency communication devices including their accessories should keep a minimum distance of 30 cm to Elipar DeepCure-L. Otherwise, the performance characteristics of the device may be reduced.
 23. Prior to each use of the device ensure that the emitted light intensity is sufficient to safely guarantee polymerization. Check the light guide and the light guide mounting hole to make sure they are clean. If necessary, the light guide mounting hole and the light guide can be cleaned as described in the section "Maintenance and Care" (see also the section "Measurement of Light Intensity").
- Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.

Place the handpiece on a flat surface when not in use. The handpiece can be connected to the charger between applications so that the battery is charged. **The handpiece must be connected to the charger** at the latest when the battery power level display glows red steadily so that the battery is charged (cf. "Battery Power Level Display on Handpiece").

The device is shipped with a light guide with 10 mm diameter. It is not permissible to use the light guides of other devices.

The handpiece is equipped with a sleep mode to minimize the device's energy consumption. The handpiece switches to sleep mode if it is not used for a period of about 5 minutes or the incorrect charging voltage is detected.

The charger uses a maximum of 0.2 W when ready for operation.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding Instructions for Use.

Product Description

Elipar DeepCure-L is a high-performance LED light source for the polymerization of light-curing dental materials. The device has two components: a wireless handpiece with a built-in battery which can be replaced by 3M customer care and a charger. The device is a medical electrical device in accordance with IEC 60601-1 and is available as a tabletop device. Wall-mounting is not possible.

Wall-mounting is not possible.

In comparison with conventional light-curing devices, Elipar DeepCure-L features excellent beam collimation and a uniform beam profile, directing more of the light energy to the restoration being polymerized and producing a deep, uniform, and complete curing of the restoration.

The light source is a high-performance light diode (LED). The beam emerging from the device covers the light wavelength range of 430 to 480 nm relevant, for instance, for camphor quinone products and is suitable for use with the majority of light-curing dental materials, including materials for fillings, liners, core build-ups, fissure sealings, temporary restorations, and cements for indirect restorations.

See the manufacturer's information for the exposure time required for the specific dental material.

Settable exposure times:

- 5, 10, 15, 20 sec
- Continuous mode (120 sec)
- Tack-cure mode (1 sec)

Intended Purpose

Intended users: educated dental professionals, e.g., general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

Fields of Application

- Polymerization of light-curing dental materials with photo initiator for the wavelength range 430–480 nm.
 - Though the majority of light-curing dental materials are responsive in this range of wavelengths, you may wish to contact the manufacturer of the material in question.

Technical Data

Charger Model CPS008050100

Operating voltage: 100–240 V 50/60 Hz

Nominal consumption: 0.2 A max

Dimensions without

country-specific adapter: Length: 65 mm

Width: 40 mm

Depth: 31 mm

Weight: 75 g

Classification: Protection class II, 

Manufacturer: Click Technology Co., Ltd.

Model: CPS 008050100

Charger Model GTM96060-0606-1.0

Operating voltage: 100–240 V 50/60 Hz

Nominal consumption: 0.6 A max

Dimensions without country-specific adapter: Length: 74 mm
Width: 44 mm
Depth: 37 mm

Weight: 107 g

Classification: Protection class II, 

Manufacturer: GlobTek, Inc.

Model: GTM96060-0606-1.0

Handpiece

Power supply: Lithium-ion battery, nominal voltage 3.7 V, capacity 2300 mAh

Utilizable wavelength range: 430–480 nm

Wavelength peak: 444–452 nm

Light intensity (between 400 and 515 nm): 1470 mW/cm² -10%/+20% (independent of battery power level)

Light emission area: 60–65 mm² (optically active)

Intermittent operation: The device has been designed solely for short-term operation. Typical operating time at room temperature (23 °C): 7 min, at 40 °C ambient temperature: 1 min on, 15 min off (cooling-off period)

Total exposure time with new, fully charged battery: Typically 120 min

Dimensions: Diameter: 28 mm
Length: 270 mm

Weight: 180 g (incl. light guide)

Charger and Handpiece

Time to charge empty battery: Approx. 2 h


Operating temperature: 10 °C up to 40 °C / 59 °F up to 104 °F

Relative humidity: 30% up to 75%

Atmospheric pressure: 700 hPa up to 1060 hPa

Transport and Storage Conditions:


Ambient temperature range: -20 °C up to +40 °C / -4 °F up to +104 °F

Relative humidity: 30% up to 75%

Atmospheric pressure: 700 hPa up to 1060 hPa

Subject to technical modification without prior notice.

Installation of the Unit

Factory Settings

The factory settings of the device are as follows:


- 10 sec exposure time


Initial Steps

Charger

- Select the plug adapter specific to the country and place it on the charger.

Charger Model CPS008050100




- ▶ Then attach the light guide to the handpiece until it noticeably locks into place (see Section “Removing and Inserting the Light Guide from/into the Hand-piece”).
- ▶ If the device malfunctions, insert the charger plug into the charging socket of the handpiece. The device will reset itself and can then be used again.

Battery Charging

- ▶ The device contains a powerful lithium-ion rechargeable battery. This type of battery does not have any memory effect and can therefore be recharged at any time by inserting the charger plug into the charging socket of the handpiece (see the section “Battery Power Level Display on Handpiece”).
- ▶ Before using the handpiece for the first time, connect it to the charger for a period of about 2 hours so that the new battery is completely charged for the first time. Make sure the colors line up correctly.

Charger Model GTM96060-0606-1.0


The green status light on the handpiece blinks while the battery is charging. The green status light glows steadily when the device is fully charged. **As a safety precaution, light-curing is not possible while the device is charging.**

- ▶ Remove the plug adapter by pushing the slider switch in the direction of the arrow.


Light Guide/Handpiece

- ▶ Place the glare shield on the front of the device.
- ▶ Autoclave the light guide prior to first use.


Battery Power Level Display on Handpiece

Status LED	Operating status	
	Handpiece without charger	Charger connected
Steady green light	Handpiece ready for operation	Charging has been completed
Flashes green	—	Battery is being charged
Steady red light	Low battery charge	Problem during charging
Flashes red	Battery fully discharged, exposure cycle will be completed or, if in continuous mode, stopped	Problem in charging, battery is defective or cannot be charged

Operation

Selection of Exposure Time

Exposure time options: 5, 10, 15, 20 sec, continuous mode (120 sec), tack-cure mode (1 sec).

- ▶ See the instructions for use for the specific dental material when selecting the exposure time.
- ▶ The indicated exposure times assume that the light guide is held at the exact position over the material being polymerized.
- ▶ If the distance between the light guide and the restoration is increased, the exposure time must be adjusted accordingly because the light intensity weakens (see graph).

Select the exposure time by pressing the TIME button.


- The selected exposure time is indicated by the 4 green LEDs.
- Each time the button is briefly pressed, the setting advances to the next (higher) value. All 4 green LEDs will be turned on for a setting of 20 sec. Pressing the button again will turn off all of the LEDs and enable the continuous mode.
- The display advances through the available settings if the button is kept depressed.
- While exposure is ongoing, the button for selection of the exposure time is inactive.

Activating and Deactivating the Light

- ▶ Briefly press the START button; the light will turn on.
 - The LEDs first indicate the set exposure time; all 4 LEDs light up for 20 sec. Every 5 sec, as the time runs down, the LEDs will turn off one at a time; at 15 sec remaining time, 3 LEDs will still be on, at 10 sec remaining time 2 LEDs, etc.
 - The LEDs do not come on at all in continuous mode; an audible signal is emitted every 10 seconds.
- ▶ If desired, the light can be turned off by pressing the START button again before the exposure time is over.
- ▶ Holding down the START button activates the tack-cure mode: the device emits a single short light pulse which enables the defined curing of Protemp™ Crown temporary restorations or a light-curing cement excess (e.g., RelyX™ Unicem) to enable easy removal.

Inserting and Removing the Light Guide from/into the Handpiece

- ▶ Attach the light guide to the handpiece by inserting it with a slight rotation until it clicks firmly into place and there are no gaps between the neck of the light guide and the handpiece (see picture).

- ▶ Remove the light guide by turning slightly and pulling it out of the handpiece.

Positioning the Light Guide


- ▶ Rotate the light guide into the desired position for polymerization.
- ▶ To make full use of the light intensity provided, place the light guide as close to the filling as possible. Avoid directly contacting the filling material.
 - Keep the light guide clean at all times to obtain full light intensity.
 - **Damaged light guides substantially reduce the light intensity and must be replaced immediately! Sharp edges may cause serious injury!**

Testing of Light Intensity

Depending on the layer thickness of the filling material being cured, it is possible to check the function of the curing light using curing disks for composites:

- ▶ Place the curing disk on a mixing pad.
- ▶ Fill your preferred composite into a bore of the curing disk at least twice as deep as the recommended layer thickness of your composite.
- ▶ Cure the composite in the curing disk for the time recommended by the composite manufacturer.
- ▶ Scrape off soft material from the bottom of the cured material with a plastic spatula.
- ▶ The solid thickness of the cured material in the curing disk **divided by two** is the layer thickness which can be properly cured.
- ▶ If you have doubts about the correct function of your curing light, clean the light guide and the protection glass in the light guide mounting hole and ensure proper fit of the light guide in the handpiece and repeat the intensity test.
- ▶ If you have still doubts about the correct function of your curing light, contact 3M service.

The curing disk must only be used to test the function of the curing light! For clinical depth of cure, please refer to the Instructions for Use of the filling material.

Commonly available devices for measuring light intensity can also be used; their measurement values should not be regarded as absolute values. When such devices are used, we recommend recording the intensity of the curing light before its first use and measuring it again at regular intervals so that any decrease in light intensity can be detected. The light intensity testing unit in an Elipar DeepCure-S base station has a device for measuring intensity on a percentage basis.

Recommended Cure Techniques with 3M Adhesives

- Rotate the light guide into the desired position for polymerization.
- A barrier sleeve may be used. Note that sleeve can impact power and irradiance of a curing light. It is critical to ensure the sleeve is powder free and the seam is not over the tip of the light guide and that the sleeve fits well. Intermediate level disinfection (manual disinfection with ready-to-use disinfectant wipes) is indicated when a sleeve is used, otherwise the light guide should undergo high level disinfection (manual disinfection using disinfectant bath, see "Reprocessing the Light Guide").

The following sleeves have been tested:

Sleeve	Impacts Power or Irradiance
Cure Sleeve® Tip Covers (short) – Kerr	No
Complete Curing Light Sleeves-Demi – Kerr	No
TIDiShield™ Custom Fit Curing Light Sleeves – TIDI Products	No

Sleeve	Impacts Power or Irradiance
SmartLite® Max Lens Cover Sleeve – Dentsply	No
Curelastic™ Cure-Light Sleeve (Steri-Shield)	Yes

Sleep Mode

The handpiece begins charging automatically when connected to the charger (green status light blinks) if the battery requires charging. If the charging voltage is not correct (e.g., there is dirt on the contacts in the charger socket or the charger), the handpiece goes into sleep mode. If the handpiece is not connected to the charger and is not used for a period of about 5 minutes, it also goes into sleep mode. In this operation mode, all of the displays and signals of the handpiece are turned off so that power consumption is reduced to a minimum. To terminate the sleep mode, press the START button.

- The sleep mode termination signal (two short audible signals) is emitted, indicating that the handpiece is ready for operation; the handpiece displays the latest selected exposure mode and time settings.

Audible Signals – Handpiece

An audible signal is emitted

- every time a button is pressed,
- every time the light is turned ON,
- 1 time after 5 sec exposure time, 2 times after 10 sec, 3 times after 15 sec.
Exception: in continuous mode, an audible signal is emitted every 10 seconds.

Two audible signals are emitted

- every time the sleep mode is terminated by pressing the START button,
- every time the light is turned OFF.

A 2-sec error signal is emitted, if

- the handpiece overheats,
- the battery lacks sufficient charge.

The audible signals from the handpiece can be turned off (except for the 2-sec error signal). Follow these instructions to turn them off. Plug the charger into a functioning power outlet, then press and hold down the TIME and START buttons at the same time with one hand. With the other hand, connect the charging cable

from the power supply to the charging socket of the handpiece. An audible signal confirms that there has been a change from “Audible signals activated” to “Audible signals deactivated”. Release the two buttons and disconnect the handpiece from the charger. Repeat the above procedure to activate the audible signals.

Troubleshooting

Error	Cause ► Solution
The status display on the handpiece glows red steadily.	The remaining battery charge is adequate only for typical 5 10-sec exposure cycles. ► Connect the handpiece to the charger and recharge the battery.
The status display of the handpiece flashes red. The ongoing exposure is interrupted (light off signal is emitted) followed by a 2 sec-error signal; the handpiece switches to “sleep” mode and resists further activation.	The battery lacks sufficient charge. ► Connect the handpiece to the charger and recharge the battery.
The status light on the handpiece blinks red whenever the handpiece is connected to the charger.	Charging problem. The battery is defective or at the end of its useful life. ► Contact 3M Service.
The handpiece has not been used for a long time and now it cannot be turned on.	There is not enough charge in the battery to turn on the handpiece. Connect the handpiece to the charger and recharge the battery.
The handpiece does not respond to the pressing of either button.	Software crash possible. ► Plug the included charger into an outlet and connect it to the handpiece. This causes the light-curing device to reset itself.

Error	Cause ► Solution
While the handpiece is connected to the charger, pressing the START button will not start the light emission.	The handpiece is connected to the charger. Light-curing is prevented as a safety precaution. ► Disconnect the charger from the handpiece and restart the light emission.
The light emission does not start when the START button is pressed; an error signal sounds for 2 sec.	The information about intermittent operation under the header "Technical Data – Handpiece" has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.
During light emission in continuous mode, an error signal sounds for 2 sec, the emission is stopped, and the handpiece changes to sleep mode.	The information about intermittent operation under the header "Technical Data – Handpiece" has not been observed. The handpiece has become overheated in the course of use. The handpiece can be used again once it has cooled down. ► Allow the handpiece to cool for 3 minutes, then start the next exposure by pressing the START button.
The light intensity is too low.	► Clean the light guide and the protecting glass in the light guide mounting hole (please refer to "Cleaning the Light Guide").

Error	Cause ► Solution
The dental material does not cure completely.	► Clean the light guide and the protecting glass in the light guide mounting hole (please refer to "Cleaning the Light Guide"). ► Check to see that the correct light guide has been attached.
The light guide cannot be attached to the handpiece.	► The light guide is not designed for use with the Elipar DeepCure-L.

Maintenance and Care

The Elipar DeepCure-L device is maintenance-free. No periodic maintenance is required. See the information contained in this chapter to secure problem-free operation.

Care of the Handpiece

- Use only the charger included with the product. The use of other chargers may damage the battery cells or result in inadequate charge.

Do not immerse the handpiece in water or incinerate. Please also observe the chapter on "Safety".

Reprocessing the Light Guide

The light guide is not sterile when delivered and must be reprocessed before being used for the first time and after each patient use.

Follow manufacturer's instructions for preparation, temperatures, contact time and rinsing off the cleaning¹ and disinfecting solution² exactly.

An enzymatic detergent should be used as the cleaning solution.

Material Resistance

Make sure that the cleaning and disinfectant agents you have chosen do not contain any of the following materials:

- Organic, mineral, and oxidizing acids (minimum acceptable pH value 5.5)
- The pH value of the cleaning and disinfecting solution should be between pH 7 and 11
- Oxidation agents (e.g., hydrogen peroxide)
- Halogens (chlorine, iodine, bromide)
- Aromatic/halogenized hydrocarbons

Observe the manufacturer's information about the cleaning and disinfecting agents.

The light guide must not be exposed to temperatures higher than 135 °C (275 °F).

The light guide has been tested for up to 500 sterilization cycles.

Limitations on Reprocessing

This medical device is safe to use if it is reprocessed according to these instructions and no visible damage is seen.

Inspect the device carefully before each use and discard if there is any damage.

Point of use / Pre-Treatment

The contaminated light guide must be placed in an approved, sealed container during transport from the point of use to the reprocessing area.

The pre-treatment must be carried out before either automatic or manual cleaning and disinfecting. Immediately after using (within a maximum of 2 hours), remove gross contaminations from the light guide.

For this, rinse the light guide thoroughly under running water (drinking water quality as a minimum) for at least 10 seconds.

Use a soft brush or a soft cloth to manually remove any remaining contamination. Remove polymerized composite with alcohol; a plastic spatula may help in removing the material. Do not use any sharp or pointed tools to protect the surface of the light guide from scratching.

Cleaning and Thermal Disinfecting (Automated Washer-Disinfectors / CDD (Cleaning and Disinfection Device))

- ▶ After the light guide has been pre-treated, place it in the wire or sieve basket of the washer-disinfector. Make sure that the devices do not touch one another in the washer-disinfector.
- ▶ Remove the light guide from the washer-disinfector unit at the end of the drying cycle and store it under low germ count conditions.
- ▶ Examine the light guide for visible impurities. If any impurities are visible, repeat the cycle in the washer-disinfector.

Automated washer-disinfectors that meet the following requirements are suitable for use:

- The device generally has a certified effectiveness in accordance with ISO 15883. The washer-disinfector is regularly serviced and tested.
- If available, use a tested cycle for thermal disinfection (AO value > 600 or, for older devices, at least 5 min at 90 °C/194 °F).

- The water used for rinsing is minimum drinking water quality. Use low-germ and low-endotoxin water for final rinsing.
- The air used for drying should be filtered (oil-free, with low microbiological and particulate contamination).

Note: The evidence of suitability for automatic cleaning and disinfection took place using the Miele G7336 CD washer-disinfector (Gütersloh) and the cleaning agent Thermosept RKN-zym (Schülke & Mayr). The program D-V-THERMO0905, without drying time was used under worst case conditions.

Cleaning steps	Parameters	
Pre-rinse	Temperature:	10 ° ± 2 °C
	Time:	1 min
Cleaning	Cleaning temperature:	45 ° ± 2 °C
	Time:	5 min
	Cleaning solution:	Thermosept RKN-zym
	Concentration:	0.3% (3 ml/l)
Final rinse	Temperature:	10 ° ± 2 °C
	Time:	2 min
Thermal disinfection	Temperature	90 ° ± 2 °C
	Time:	5 min

Cleaning and disinfection of the light guide by wiping is only permitted in combination with a corresponding protective sleeve

Manual Cleaning

- ▶ Clean the light guide after every application for 30 seconds with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds (QAV)) or equivalent cleaning agents.
- ▶ Follow the cleaning agent manufacturer's directions for use.
- ▶ Clean with the wipes until no contamination is visible.
- ▶ Then rinse with water (drinking water quality as a minimum) and dry with a soft, lint-free cloth.

Manual Disinfection (Intermediate Level Disinfection)

- ▶ An intermediate level disinfection is possible only if the light guide is covered with a protective sleeve during use.

- ▶ Disinfect the light guide with ready-to-use disinfectant wipes (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds (QAV)) or equivalent disinfection agents. Follow the disinfection agent manufacturer's directions for use, observing in particular the contact times.
- ▶ Then rinse with water (drinking water quality as a minimum) and dry with a soft, lint-free cloth.

Manual Cleaning and Disinfection the Light Guide by Immersion Bath

Manual Cleaning

- ▶ A neutral enzymatic cleaning agent is recommended (e.g., Cidezyme/Enzol from Johnson & Johnson).
- ▶ Place the light guide for the specified application time into the solution, making sure that it is completely covered.
- ▶ Cleaning takes place in the immersion bath with a soft brush for at least 5 minutes.
- ▶ Then remove the light guide from the solution, thoroughly rinse with low-germ-count water (drinking water quality as a minimum), and dry with a lint-free cloth.

Manual Disinfection (High Level Disinfection)

- ▶ A disinfecting agent containing o-Phthalaldehyde (e.g., Cidex OPA by Johnson & Johnson) is recommended. It is compatible with the light-curing device.
- ▶ To disinfect, place the cleaned light guide for the specified application time (≥ 12 min) into the solution, making sure that it is completely covered.
- ▶ The instruments in the disinfection bath must not touch each other during this process.
- ▶ Then remove the light guide from the solution, thoroughly rinse with low-germ-count water (drinking water quality as a minimum) for at least 3 minutes, and dry with a lint-free cloth.
- ▶ Then check the light guide (see "Check").

¹ Cleaning agent Thermosept RKN-zym (Schülke & Mayr) for the washer-disinfector and Cidezyme® / Enzol® and CaviWipes™ were used for manual cleaning.

² Cidex OPA, ready-for-use solution, was used for validation of the high level disinfection and CaviWipes™ for validation of intermediate level disinfection.

Packaging for Sterilization

Package the light guide in single wrapped sterilization pouches (one-way packaging or double packaging).

Sterilization

Effective cleaning and disinfection are absolutely essential requirements for effective sterilization.

A fractionated vacuum procedure (at least 3 vacuum stages) or a gravity procedure (with adequate product drying) must be used for the sterilization. The use of the less effective gravity procedure is permissible solely if the fractionated vacuum procedure is unavailable. The user is solely responsible for selection of the parameters required for the gravity procedure; the user must therefore calculate or validate these parameters in accordance with the loading configuration. The drying time actually required is also directly dependent on parameters that are in the sole responsibility of the user (loading configuration and density, sterilizer condition, etc.) and must consequently be determined by the user.

- Steam sterilizer in conformity with DIN EN 13060/DIN EN 285 or ANSI AAMI ST79 (for USA: FDA clearance)
- Maximum sterilization temperature 135 °C (275 °F)
- The following cycles were validated for the light guide:

	Fractional prevacuum		Gravity displacement
Temperature	132 °C (270 °F)	134 °C (273 °F)	121 °C (250 °F)
Sterilization time	4 min	3 min	20 min
Drying time	30 min	30 min	30 min

Check

Before using the light guide again, check it for damaged surfaces, discoloration, and contamination; do not use damaged light guides. If the light guide is still contaminated, repeat all reprocessing steps.

Clean and Disinfect Handpiece and Glare Shield

- Cleaning and disinfecting solutions must not enter the unit.
- Wipe the surfaces with ready-to-use disinfectant wipes³ (e.g., CaviWipes™ based on the active substances alcohols and quaternary compounds (QAV) or equivalent disinfection agents). Subsequently inspect the surfaces for visible impurities. If any impurities are visible, repeat the cleaning procedure.
- Solvents or abrasive cleaners may not be used as they can damage the components.
- Wipe the device surfaces with a tuberculocidal disinfectant wipe³ for disinfection and keep it wet according to the disinfection agent manufacturer's directions for use, particularly observing the contact time.
- Disinfectants must not be sprayed directly onto the device surface.
- Wipe the device surfaces with a soft cloth moistened with water (drinking water quality as a minimum) and

dry with a soft, lint-free cloth. If you do not rinse well with water, the components will be damaged.

Additional notes: Ensure that the charge contact pins stay dry and are not touched by any metallic or greasy parts. Do not bend the charge contact pins while they dry. Wet charge contact pins cause an operating error (malfunction message: the LED showing the operating status of the charger flashes red).

³ CaviWipes™ were used for validation of the cleaning and disinfection procedure.

Clean Handpiece and Glare Shield

Clean all components with a soft cloth and, if necessary, a mild cleaning agent (e.g., dish-washing detergent). Solvents or abrasive cleaners can damage the components.

Cleaning agents must not enter the device.

- ▶ To disinfect all components, spray the disinfectant on a towel and use it to disinfect the device. Do not spray the disinfectant directly on the handpiece.
 - Disinfection agents must not enter the device!
- ▶ Dry residual disinfectants with a soft and fluff-free cloth as they damage the plastic components.
- ▶ Make sure that disinfectants do not come into contact with the charging socket on the handpiece because this could impair proper charging.

If necessary, ask the manufacturer of the disinfectant if its constant use will damage plastic surfaces.

Clean the protection glass with a soft and fluff-free cloth. Beware of scratches!

Storage of the Handpiece during Extended Periods of Non-Use

- ▶ If the handpiece will not be used for a number of weeks – e.g., during vacation – charge the battery beforehand or connect the handpiece to the charger for this time. A safety switch within the battery prevents a total discharge.

Discharged or nearly discharged batteries must be recharged as soon as possible.

Return of Old Electric and Electronic Equipment for Disposal

1. Collection

Users of electric and electronic equipment are required to collect their old equipment separately from other waste in accordance with the regulations of the specific

country. Old electric and electronic equipment must not be disposed of with unsorted household waste. This separate collection is a prerequisite for recycling and reprocessing as an important method for preserving environmental resources.

2. Return and Collection Systems

When your Elipar DeepCure-L is no longer usable, do not dispose of the device with household waste. 3M Deutschland GmbH has set up special disposal facilities to handle the equipment. Details about the procedure for the specific country can be obtained from the pertinent 3M subsidiary.

3. Removing the Battery for Disposal

To dispose of the battery, disconnect the Elipar LED curing light from the charger, remove the two screws on each side of the charging socket and push the lower half of the housing towards the back, away from the upper half of the housing. Use a suitable tool to cut the connecting wires between the battery and the circuit board and remove the battery for disposal in accordance with 1. and 2.

Customer Information

No person is authorized to provide any information that deviates from the information provided in this instruction sheet.


Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Symbol Glossary

Reference Number and Symbol Title	Symbol	Description of Symbol
ISO 15223-1 5.1.1 Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.
ISO 15223-1 5.1.3 Date of Manufacture		Indicates the date when the medical device was manufactured.
ISO 15223-1 5.1.5 Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified.
ISO 15223-1 5.1.6 Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified.
ISO 15223-1 5.1.7 Serial number		Indicates the manufacturer's serial number so that a specific medical device can be identified.
ISO 15223-1 5.3.7 Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed.
ISO 15223-1 5.3.8 Humidity limitation		Indicates the range of humidity to which the medical device can be safely exposed.
ISO 15223-1 5.3.9 Atmospheric pressure limitation		Indicates the range of atmospheric pressure to which the medical device can be safely exposed.
ISO 15223-1 5.4.4 Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself.
CE Mark		Indicates conformity to European Union Medical Device Regulation or Directive.
Medical Device		Indicates the item is a medical device.
Follow instructions for use		To signify that the instructions for use must be followed.
Rx Only		Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1).

Reference Number and Symbol Title	Symbol	Description of Symbol
Type B applied part		To identify a type B applied part complying with IEC 60601-1.
Class II equipment		To identify equipment meeting the safety requirements for Class II equipment per IEC 61140.
Use Indoors		Indicates medical device be used indoors
Regulatory Compliance Mark		Indicates product complies with applicable Australian regulatory requirements.
Battery		Battery power level
PAP 20/21		Indicates paper material is recyclable.
Recycle electronic equipment		DO NOT throw this unit into a municipal trash bin when this unit has reached the end of its lifetime. Please recycle.
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.

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Informasjonsstatus: januar 2020

Guidance and manufacturer's declaration – electromagnetic emission		
The Elipar DeepCure-L is intended for use in the electromagnetic environment specified below. The customer or the user of the Elipar DeepCure-L should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The Elipar DeepCure-L uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. The Elipar DeepCure-L is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies	

Immunity test	Compliance to IEC 60601 test levels	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radiated RF EM fields IEC 61000-4-3	3 V/m 80 MHz - 2,7 GHz 80 % AM at 1 kHz	
Proximity fields from RF wireless communications equipment IEC 61000-4-3	see following table: Test specifications for enclosure port immunity to RF wireless communications equipment	
Electrical fast transients / bursts IEC 61000-4-4	± 2 kV for power supply lines 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Conducted disturbances induced by RF fields IEC 61000-4-6	3 V _{rms} 0,15 MHz - 80 MHz 6 V in ISM bands between 0,15 MHz and 80 MHz 80 % AM at 1 kHz	
Surges Line-to-line IEC 61000-4-5	$\pm 0,5$ kV, ± 1 kV	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, and short interruptions on power supply input lines IEC 61000-4-11	0 % U _T ; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315° 0 % U _T ; 1 cycle and 70 % U _T ; 25130 cycles Single phase: at 0° 0 % U _T ; 2501300 cycle	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Elipar DeepCure-L Dental Light Unit requires continued operation during power mains interruptions, it is recommended that the Elipar DeepCure-L Dental Light Unit be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U _T is the a. c. mains voltage prior to application of the test level..		

Test specifications for enclosure port immunity to RF wireless communications equipment

Test frequency (MHz)	Band (MHz)	Service	Modulation	Maximum power (W)	Distance (m)	Immunity test level (V/m)
385	380 – 390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460 FM	± 5 kHz deviation 1 kHz sine	2	0,3	28
710 745 780	704 – 787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
810 870 930	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
1 720 1 845 1 970	1 700 – 1 990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
2 450	2 400 – 2 570	Bluetooth, WLAN, 802 .11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5 240 5 500 5 785	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9

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Health Care Business

Carl-Schurz-Str. 1

41453 Neuss - Germany

3M Customer Care/MSDS Information: U.S.A. 1-800-634-2249.

3M, ESPE, Elipar, Protemp and RelyX are trademarks of 3M or 3M Deutschland GmbH.
Used under license in Canada. © 2020, 3M. All rights reserved.

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3M ou 3M Deutschland GmbH.