Gebrauchsanweisung Instructions for use Bio-Ink®-Flow Sore Spot Indicator - Set BK 205

Hersteller | Manufacturer

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ENGLISH

General description and purpose of the medical device: Green color paste for marking sore spots on dentures and partial dentures. For professional use only.

Indication

Bio-Ink® Flow is indicated for marking sore spots on the gum and subsequent marking on the denture. Sore Spots also known as decubitus can occur frequently with new dentures or partial dentures. Also a degeneration of the alveolar ridge can cause sore spots as well as unbalanced occlusal forces on the denture. The adaption of the oral mucosa to the new denture often causes inflammatory lesions. Bio-Ink® Flow is a highly viscous paste. By punctiform, laminar or linear application to the mucosa the Palatal Vibrating Line and the mucobuccal fold can be marked. If necessary, the denture can be adjusted.

Contraindication

A direct application onto highly inflamed areas or purulent lesions of the oral mucosa areas should be avoided.

Side Effects

Possible side effects could be allergic reactions.

Notes

- The color paste can be applied directly with the single use cannulas and spread with the single use applicator (e.g. Magic-Brush™) to avoid cross-contamination.
- Alternatively, the paste can be used with a single use applicator (e.g. Magic-Brush™) and then applied to the surface to be tested.
- The mucous membranes, palate and mandible must be dried with clean pulp or cotton roll, previously.
- The set contains components (Magic-Brush™ applicator, single use cannulas and a Multi-Hold™ holder for brush). Except the Multi-Hold™ holder for brush, all components are for single use only and must not be re-processed.
- The product must be stored in the sealed syringe at temperatures between 0°C and 40°C.

Package content

- 1 x Bio-Ink® Flow Paste, Syringe 1ml
- 10 single use Applicator Needles (cannulas)
- 1 holder for brush (Multi-Hold™ for Magic-Brush™)
- 50 brushes Magic-Brush™ (single use product)
- Instructions for use

Preparation

- Before use, the mucosa membrane, the palate and the mandible must be dried carefully with a pulp or a cotton roll so that the color paste does not run.
- If necessary, clean the holder for the applicator.
- If necessary, fit the applicator into the holder.

Application

Before use make sure that the mucosa, palate and alveolar ridge are dried carefully with a cellulose pulp or a cotton roll so that the color paste does not run. There should be no saliva on the surface to be tested. The denture must also be cleaned and dried carefully beforehand. Insert the applicator (Magic-Brush™) into the provided holder (Multi-Hold™).

The color paste can be applied directly with the disposable cannulas and spread with the disposable applicator (e.g. Magic-Brush™). Alternatively, the color can be used with the disposable applicator (z.B. Magic-Brush™) and then applied to the surface to be tested.

To facilitate handling, the applicator can be angled at the tip.

Effect/Diagnosis

The color marking of the sore spot will be transferred to the denture when the denture is seated into position. When the denture is removed the spot is visible on the surface of the denture. If necessary the denture can be adjusted by the dentist.

Cleaning (of the surface to be tested)

After control, the color marking can be easily removed with water and a suitable surfactant (dishwashing detergent).

Reprocessing

In principle, the following should be observed: „Requirements for hygiene in the processing of medical devices Recommendation of the Commission for hospital hygiene and infection prevention (KRINKO) at the Robert Koch Institute (RKI) and the Federal Institute for Drugs and Medical Devices (BfArM)“, Bundesgesundheitsblatt 2012-55: 1244-1310.

- The instrument (holder) must be reprocessed immediately after each treatment (within a maximum of 2 h).
- Thoroughly clean the holder by hand with a small brush with firm bristles.
- Note: The cleansing should be done in a water bowl without further additives below the water surface in order to achieve adequate cleansing of the instrument parts whilst avoiding protein fixation as well as to protect the environment from contamination due to splash water.
- Immediately rinse the holder with water (minimum drinking water quality).
- Place the instrument in a standard cleaning and disinfection bath. Examples: Becht Premium Konzentrat, Dürr Dental ID 213 Instrument disinfection, Pluradent instrument bath PLULINE, Schülke & Mayr gigasept® instru AF (see also „List of disinfectants and disinfection methods tested and approved by the Robert Koch Institute“ or disinfectant list VAH).
- Note: Formaldehyde containing cleaning and disinfecting agents may only be used after sufficient cleaning in order to avoid protein fixation.
- Note: The instructions for use of the manufacturer of the cleaning and disinfecting agent must be strictly observed. In particular, the concentrations to be used and the reaction times must be observed!
- Note: For machine cleaning, the manufacturer's instructions for use must be followed carefully!
- Final rinse of the instrument parts with water (minimum drinking water quality, recommended: demineralized water with microbiological quality according to drinking water).
- Dry.
- Visual inspection for damaged surfaces, chipping and impurities. Damaged instruments must be discarded (For a limited restriction of reuse, see the chapter „Reusability“). If residual contamination still exists, the entire cleaning process must be repeated with all steps (Cleaning, intermediate rinsing, disinfection, final rinsing and drying).
- The instrument must be free of any residue and dry before further processing.
- Maintenance of the instruments is not required.
- The instrument is packaged and sealed in a sufficiently large single sterilization bag according to EN 868-2 ff ISO 11607 (suitable for steam sterilization). Observe the instructions of the manufacturers of the sterilization bags and the sealing machines and the current normative requirements.
- Sterilization must be carried out in a validated process with moist heat in an autoclave according to DIN

EN 13060 Type B or DIN EN 285 resp. ANSI AAMI ST79. Observe the manufacturers operating instructions of the autoclave.

- Sterilize the instrument with humid heat (saturated water vapor) using a pre-vacuum procedure for 5 minutes at 134°C.
- The recommended storage life for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling.

Reusability

- Frequent reprocessing restricts the life of the plastic instrument. The instrument must be disposed of in case of cracks in the plastic.
- It is not permitted to reuse single use items (applicator (Magic-Brush™) as well as single use cannulas). Safe use cannot be guaranteed if these products are used again, since there is an infection risk and / or the safety of the products is not liable.
- The use of damaged or soiled instruments is the responsibility of the user.
- In case of non-compliance, all liability is excluded.

Disposal Instructions

Bio-Ink®-Flow consists of a green food color and natural waxes, which are approved for food, cosmetics and pharmaceutical applications. All components of the kit can be disposed with the usual practice waste.

Notification of Incidents

Serious adverse events occurring in connection with the product must be reported to the manufacturer and the competent authority of the Member State.

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