DENTA-FAST NF

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Directions for use

Indications:
For repairs, extending dentures and reattaching teeth, relinings.

Advantages:
Short processing time, the acrylic dough can be processed immediately.

Mixing ratio:
Powder and liquid at a ratio of 5 : 3
• 5 parts by weight powder
• 3 parts by weight liquid

Processing:
• Apply alginate separating agent to the plaster.

• Add the powder to the liquid and spatulate the Denta-Fast NF for 15 seconds to mix it homogeneously.

• To mix the required amount of Denta-Fast NF, pour liquid into a mixing vessel and sprinkle powder in until it is just below the surface of the liquid. Mix for approximately 15 seconds.

• The acrylic dough has a runny consistency, form drops and can be processed immediately.

• When reattaching teeth, as is usual when using cold-curing resins, they must be roughened, provided with undercuts and wetted with monomer or adhesive agent.

• Once the surface of the denture has been roughened, the runny dough should be dripped onto it and spread evenly. Once the resin has set, the denture can be ground and polished immediately.

• Once the dough has attained a plastic consistency, it must not be touched or contoured.

• Larger amounts of Denta-Fast NF should be applied in layers to prevent porosity. The only method of ensuring reliable, non-porous setting is to cure the material under pressure.

Polymerization:
Pressure vessel for 15 minutes at a pressure of 2-5 bars and 45°C water temperature.

Storage:
At 10-25°C, please note the expiry date!

Important:
Should the room temperature deviate from 21°C: Higher temperatures curtail the processing time and lower temperatures prolong it.

Packaging:
Powder: 500g or 1000g
Liquid: 500ml or 1000ml

Danger:

Avoid prolonged skin contact with non-polymerized material an inhalation of monomer vapours.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful the side-effects has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend no to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Troubleshooting

Problem	Cause	Remedy
Acrylic does not fully cure	Incorrect mixing ratio	• Adhere to the instructions for use • Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	• Incorrect or poor application of separating agent • Model not soaked	• Apply ample separating agent • Check the separating agent if necessary • Soak the model for approx. 10 min. in handhot water
Acrylic surface white / porous	Recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Color differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic / denture has whitish streaks	Incorrectly mixed / too much powder used	Thoroughly spatulate the acrylic dough / adhere to the mixing ratio

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/ Germany · Tel.: +49 -6003-814-0 · Fax: +49 -6003-814-901
Vertrieb durch Müller & Weyandt GmbH · Reichardsweide 40 · D-63654 Büdingen · Tel.: +49 (0) 60 42 - 88 00 88