DENTA-PRESS

Stand:250521

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1. Allgemeine Produktbeschreibung

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Processing instructions

1. General product description

- Cold-curing methyl methacrylate-based resin.

2. Composition

- Powder: Bead polymer of polymethyl methacrylate (PMMA) Pigments, initiators
- Liquid: Methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilizers
- M+W Denta-Press is cadmium-free.

3. Indications

- For the production of partial and full dentures using the flow technique
- Production of the acrylic portion of implant-borne dental appliances and model casting dentures using the flow technique
- Repairs, extensions and relinings of existing restorations using the flow technique
- For the production of splints
 - o Mouthguards
 - o Bite splints

4. Contraindications / Interactions with other agents

- Do not use in the case of a known allergy to one of the components.

5. Available colours

- The powder is available in various colours, in some cases with veins.

6. Application / Preparation

a. Preparatory work

Produce the plaster cast and wax-up of the teeth as usual.

The fabricated tooth set in wax is secured with plaster or silicone (matrix/flask). Then the wax is eliminated as usual.

To prevent artificial teeth from breaking off, they must be roughened, have undercuts made, and be wetted with monomer or Connector.

The surfaces of the plaster models are separated with an alginate separating agent (Separator) according to manufacturer's instructions.

b. Mixing ratio

10 parts by weight powder, 7 parts by weight liquid
Veins can be included, if desired. Please observe separate manufacturer's instructions.

c. Production of partial and full dentures as well as implant-borne appliances and model casting completions using the flow techniques

Mix powder and liquid (according to the mixing ratio) and after allowing it to thicken for approx. 15 seconds, stir with a spatula until a homogeneous mixture is formed.

Fill acrylic dough into the mould as usual.

After mixing, M+W Denta Press is pourable for approx. 2-3 minutes at a room temperature of 21°C; after 5-6 minutes, it transitions into the plastic phase.

After the start of the plastic phase, the acrylic no longer flows out of the matrix and it can be modelled.

The dental appliance must be placed in the pressure vessel after no more than 7-8 minutes.

The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45°C.

d. Repairs and relinings of existing acrylic restorations using the flow technique

Create a model or a matrix, if necessary.
Clean and etch the acrylic surfaces to be repaired,

smooth the edges and wet with monomer or Connector.

Mix powder and liquid (according to the mixing ratio) and after allowing it to thicken for approx. 15 seconds, stir with a spatula until a homogeneous mixture is formed.

Allow acrylic to run off.

The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45°C.

e. Finishing the material

After polymerisation, the dental appliance is lifted from the model and finished using suitable instruments.

Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide cutters and diamond instruments are suitable for finishing.

Sandpaper or pumice stone is used for pre-polishing.

Polishing

The material is polished with goat hair brushes, polishing paste and soft wool buffing wheels. Careful surface finishing and polishing is necessary for an optimal result.

Production of splints

a. Preparatory work

Production of a plaster model Blocking out of undercuts with suitable materials Model the splint in wax as desired. Make impression, create flow channels, and remove wax.

Separate the plaster models with suitable materials (e.g. Separator)

b. Processing

Mix powder and liquid (according to the mixing ratio) and after allowing it to thicken for approx. 15 seconds, stir with a spatula until a homogeneous mixture is formed.

Fill acrylic into the mould as usual.

These work steps must be completed after no more than 7-8 minutes.

The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45°C.

c. Finishing the material

After polymerisation, the splint is lifted from the model and finished using suitable instruments.

Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide cutters and diamond instruments are suitable for finishing. Sandpaper or pumice stone is used for pre-polishing.

Polishing

The material is polished with goat hair brushes, polishing paste and soft wool buffing wheels. Careful surface finishing and polishing is necessary for an optimal result.

d. Repairs

Create a model or a matrix, if necessary.

Clean and etch the acrylic surfaces to be repaired, smooth the edges and wet with monomer or Connector. Mix acrylic according to mixing ratio and allow it to run off. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2- 5 bar and a water temperature of 45°C.

7. Troubleshooting see table

8. Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own

experience and should only be considered guidelines. Our products are subject to continued development. We therefore reserve the right to make changes.

9. General information

In the event of deviations from 21°C room temperature:

Higher temperatures shorten, lower temperatures prolong the above times.

To further reduce the risk of intolerance reactions, we recommend placing the denture in lukewarm water for at least 12 hours prior to insertion.

10. Storage

- Storage temperature 10-25°C / 50-77°F. Close container tightly.

11. Shelf life

The maximum shelf life is printed on the label of each package. Do not use after the expiry date.

12. Adverse effects

With proper preparation and use of this medical device, adverse effects are extremely rare. Immune

reactions (such as allergies) or local discomfort, however, cannot be ruled out completely. If you become aware of any adverse effects – even in cases of doubt – please notify us.

13. Contraindications/Interactions

If a patient has hypersensitivities to a component of this product, it should not be used or used only under the strict supervision of the attending physician/dentist. The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

14. Note

Please pass all of the above information on to the attending physician if you are processing this medical device for a custom-made product. When processing this medical device, observe the applicable safety data sheets.

15. Disposal

Powder can be disposed of in household waste. Liquid should be disposed of as hazardous waste.

7. Troubleshooting / FAQ List

Problem	Cause	Remedy
Acrylic does not cure properly	Incorrect mixing ratio	– Adhere to the instructions for use – Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	– Inadequate separation – Model not soaked	– Apply sufficient amount of separating agent – Check the separating agent, if necessary – Soak the model for approx. 10 min. in lukewarm water
Acrylic forms bubbles on the fitting surface	– Model not soaked	Soak the model for approx. 10 min. in lukewarm water
Acrylic surface white/porous	Acrylic injected too late or placed in pressure vessel too late	Observe the processing times
Colour differences in the case of repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Incorrectly separated	Check the separating agent
Acrylic/denture has whitish streaks	– Incorrectly mixed – Too much powder used	– Mix the acrylic dough thoroughly with a spatula – Adhere to the mixing ratio
Acrylic sets too fast	Excessively high temperatures/ excessively long work steps	Observe the temperatures/ processing times
Splint breaks	Contouring too thin/unsuitable	Ensure sufficient layer thickness (at least 2-3 mm)

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany · Tel.: +49 (0) 6003-814-0 · Fax: +49 (0) 6003-814-901