Repair and correction possibilities with ProBase Cold

- Corrections and repairs of ProBase Hot, ProBase Cold and SR Ivocap® may be carried out with ProBase Cold by using the pouring technique. Thoroughly roughen the corresponding surfaces and wet with monomer.
- The residual monomer content after polymerizing the material according to the method described is < 4,5%.

Warnings

- The monomer contains methyl methacrylate (MMA).
- Methyl methacrylate is easily flammable and irritating (flash point +10 °C / 50 °F).
- Irritating to eyes, skin, and respiratory system.
- May cause sensitization by skin contact.
- Avoid contact of the skin with monomer or uncured material. Commercial medical gloves do not provide protection against the sensitizing effect of methacrylates.
- Do not breathe vapour.
- Keep away from sources of ignition – no smoking.
- Do not empty into drains.
- Take precautionary measures against static discharges.

Storage

- Store material in a cool, dark, well-ventilated place. Storage temperature: 2–28 °C (36–82 °F).
- Do not use the materials after the indicated date of expiration.
- Keep out of the reach of children.

Date information prepared: 08/2012

The material has been developed solely for use in dentistry. Processing should be carried out strictly according to the Instructions for Use. Liability cannot be accepted for damages resulting from failure to observe the Instructions or the stipulated area of application. The user is responsible for testing the material for its suitability and use for any purpose not explicitly stated in the Instructions. Descriptions and data constitute no warranty of attributes and are not binding.

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Pulver:

Polymethylmethacrylat, Weichmacher, Benzoylperoxid, Katalysator, Pigmente

Flüssigkeit:

Methylmethacrylat, Dimethacrylat, Katalysator

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