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☞ Instructions for use

Applications:

- Block-out material to balance under-
lying dry plaster stumps and undercuts
- Block-out material to fabricate a bleach-
ing splint

Use in fabricating bleaching splints:

1. Apply M+W Bloc Out once the model is
completely dry. It is applied to the labial sur-
faces of the tooth to be bleached and subse-
quently serves as a reservoir for the gel. Gen-
erally these are teeth 15 – 25 and 35 – 45.
Of course, fewer, more or individual teeth
may be prepared.

Important: Apply the M+W Bloc Out ma-
terial approx. 0.5 mm thick and keep a
distance of 0.5 – 1 mm from the gingival
margin. Do not cover occlusal surfaces with
the material!

2. After curing, thoroughly remove the
inhibition layer on the surface. A cloth
dampened with alcohol serves this pur-
pose. Do not use acetone or similar or
similar chemicals.

3. Deep draw the plaster model prepared
this way using a suitable deep drawing foil
(e.g. M+W KopsiSoft bleach).

Polymerisation time:

- Commercially available halogen
polymerisation lamp (approx. 20 sec.
per tooth)
- Plasma lamp (approx. 3 – 6 sec.,
depending on the manufacturer)
- Light oven (approx. 2 – 3 min., depend-
ing on the manufacturer)

Storage:

The material should be stored at
10°C - 25°C.
Tightly close the syringe immediately
after use.

Shelf life:

The maximum shelf life is printed on the la-
bel of the respective syringe. No longer use
after the expiration date.

Warning:

This product contains substances which
may cause allergic reactions on skin
contact in sensitive persons. Reduce the
contact with these materials, especially
with the uncured material in order to
reduce the risk of allergic reactions. In
case of inadvertent skin contact, immedi-
ately wash off with soap and water. The
use of protective gloves and a non-con-
tacting technique for working is recom-
mended. Acrylates may penetrate com-
mercially available protective gloves.
In case of contact with the protective
gloves, remove and destroy them. Immedi-
ately wash hands with soap and water and
put on new protective gloves.

Warranty:

Our application information, irrespective of
whether it is issued orally, in written form or
by practical instructions, is based on our own
experience and can therefore only be consid-
ered as guidelines.

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