DENTA-IN-PRESS

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Directions for use

Cold-curing methyl methacrylate resin.

If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

Indications:

Fabricating full and partial dentures using the injection moulding technique.

Properties:

- Permanency of color
- High mechanical values
- Simple processing

Advantages gained by using M+W Denta-In-Press in combination with injection moulding system

- Optimum reproduction of the occlusion
- Precisely fitting full and partial dentures
- Can be poured – does not come into contact with the user's hands
- Rational to work with

Mixing ratio:

The powder and liquid should be mixed in the following ratios:
2 : 1 by weight

An average sized full denture requires the following amounts of material:
45cm³ powder: 17 ml liquid (volume) or 30 g powder: 15 g liquid (weight)

Processing:

- Apply alginate mould sealer to the plaster.
- After adding the powder to the liquid and allowing it to soak in for 15 seconds, spatulate the M+W Denta-In-Press thoroughly in a ceramic vessel and pour it into the cylinder, ensuring that no air bubbles are entrapped.

- Once the resin has become viscous, inject it into the flask.

- These stages must be completed within 8 minutes.

- To avoid creating porosities, the plaster surfaces in the flask must have cooled to at least 30 °C to 40 °C.

- As is usual when using cold-curing resins, acrylic teeth must be roughened, provided with undercuts and moistened with monomer or adhesive agent to prevent them becoming detached from the denture base.

Important:

Should the room temperature deviate from 23 °C (71°F): Higher temperatures reduce the working times and lower temperatures prolong them.

Polymerization:

Polymerization is carried out using the disk / spring system, with predefined pressure. At least 20 minutes must be allowed to elapse prior to devesting the M+W Denta-In-Press dentures.

Repairs - Please note:

Repairs may be carried out with Denta-In-Press, M+W Denta-Fast-NF and M+W Denta-Press. The fractured surfaces must be roughened, cleaned and moistened with monomer prior to applying the resin dough. When using a pressure/polymerization vessel M+W Denta-In-Press must be polymerized for 15 minutes at a pressure of 2 - 4 bars. Water temperature: 45 °C (113 °F).

Storage:

At 10 - 25 °C. Please observe the shelf-life.

Side-effects:

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, howe-

er, cannot be ruled out completely.

Should you learn about unwanted sideeffects – even if it is doubtful that the sideeffect has been caused by our product – please kindly contact us.

Contra-indications / interactions:

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note:

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

Caution:

Prolonged skin contact with unpolymerized material and inhalation of monomer vapours are to be avoided.

Disposal:

Dispose of powder with normal household waste. Liquids are classified as hazardous waste and must be disposed of according to local laws.

Troubleshooting

Problem	Cause	Remedy
Acrylic does not fully cure	Incorrect mixing ratio	- Adhere to the instructions for use - Adhere to the mixing ratio
Acrylic is whitish on the fitting surface	- incorrect or poor application of separating agent - model not soaked	- Apply ample separating agent - Check the separating agent if necessary - Soak the model for approx. 10 min. in luke-warm water
Acrylic form bubbles on the fitting surface	Model not soaked	- Soak the model for approx. 10 min. in luke-warm water
Acrylic surface white / porous	Recommended time for injecting acrylic exceeded or delay in inserting into the pressure vessel	Adhere to the processing times
Color differences with repairs	Incorrect mixing ratio	Adhere to the mixing ratio
Acrylic breaks when lifted	Poorly separated	Check the separating agent
Acrylic / denture has whitish streaks	Incorrectly mixed / too much powder used	Thoroughly spatulate the acrylic dough / adhere to the mixing ratio
Acrylic sets too fast	Temperatures too high / Working stages too long	Observe temperatures / working times Storage: 10 °C - 25 °C

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany · Tel.: +49 - 6003-814-0 · Fax: +49 -6003-814-901