3M Penta™ Elastomer

- Elastomer Syringe**
- Elastomer-Spritze**
- Siringue à élastomère**
- Siringa per elastomero**
- Seringa Elastomera**

Instructions for Use
Gebruachsinformation
Mode d'emploi

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- ENGLISH**

Product Description

The Penta Elastomer syringe is used for intraoral application of light-bodied and medium-bodied consistency impression materials, especially for material dispensed from the Pentamix™ mixer. The translucent syringe barrel allows the operator to determine the amount of material to be dispensed from the mixing device. Syringe barrels can be reused several times. For replacement, spare syringe barrels are available.

These Instructions for Use should be kept for the duration of product use. The product may only be used when the product labeling is clearly readable. For details on all additionally mentioned products please refer to the corresponding instructions for Use.

Intended Purpose

Intended users: educated dental professionals, e.g., general dentists, dental assistants, dental hygienists, who have theoretical and practical knowledge on usage of dental products.

Indications

- Application of light-bodied and medium-bodied consistency impression materials

Contraindications

None

Application

- Remove the syringe plunger from the syringe and attach the Penta Elastomer syringe directly to the mixing tip of the Pentamix unit for filling.
 - The translucent syringe barrel allows the operator to determine the amount of material to be dispensed from the mixing device.
- During or immediately after application dispense all material from the syringe and allow the impression material to set.

Cleaning and Disinfection of the Penta Elastomer Syringe

Caution
Observe the following requirements when selecting the cleaning and disinfectant agents that will be used:

- They should not contain the following substances: Oxidation agents (e.g., hydrogen peroxide) and oils.
- Follow manufacturer’s instructions for concentrations, temperatures, development times, and rinsing off the cleaning and disinfectant agents exactly. Use solely solutions that have been freshly prepared with water of drinking water quality. Use only a soft, clean, and lint-free cloth for drying.
- Never use metal brushes or steel wool for cleaning.
- Use only validated methods for cleaning and disinfection.
- Use solely cleaning and disinfectant agents tested for efficacy and compatibility with the medical device (e.g., FDA certification). Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals.

The plastic tube is only pre-cleaned and the outer surface is disinfected by wiping. The tube must not be immersed in a cleaning and disinfection bath.

Limitations on Reprocessing

The Penta Elastomer syringe is treated carefully and is neither damaged nor contaminated. It may be used as long as the plunger can be pushed in evenly. The user bears sole responsibility for any use of Penta Elastomer syringes that do not meet this criterion or that are damaged and/or contaminated.

General

The Penta Elastomer syringe is non-sterile when delivered and must be reprocessed according to instructions before initial use and prior to every further use.

Preparation for Decontamination

The medical device should be reprocessed as soon as possible after every use.

Reprocessing

The Penta Elastomer syringe must first be cleaned manually and, after cleaning, be manually disinfected in accordance with the following instructions.

Preparation

Disassemble the syringe into the individual parts.

Unscrew the knurled nut and application tip and discard the application tip. Use the plunger to push the set impression material out of the syringe body and clean the syringe with a clean, dry brush.

Pre-cleaning the Tube

Immediately after use (within a maximum of 2 hours), remove gross impurities. Use a clean, soft, dry brush to manually remove all visible impurities.

Pre-cleaning the Knurled Nut, Plunger, and Handle of the Penta Elastomer Syringe

Immediately after use (within a maximum of 2 hours), remove gross impurities from the knurled nut, plunger, and handle of the Penta Elastomer syringe.

- Use a clean, soft, dry brush to manually remove all visible impurities.
- The knurled nut, plunger, and handle of the Penta Elastomer syringe must be kept under running water (use potable water or drinking water quality, temperature < 35° C/95° F) and rinsed for at least 1 minute.

Manual Cleaning and Disinfection of the Tube by Wiping

The plastic tube of the Penta Elastomer syringe must not be immersed in a cleaning and disinfection bath.

Manual Cleaning of the Outer Surface by Wiping

- Clean the Penta Elastomer syringe before its first use and after every application with ready-to-use disinfectant wipes (for example CavWipes™ alcohol-quaternary ammonium disinfectant) or equivalent cleaning agents.
- Follow the cleaning agent manufacturer’s directions for use.
- Continue wiping the syringe with the wipes until no contamination is visible.

Manual Disinfection of the Outer Surface by Wiping (Intermediate Level Disinfection)

- Disinfect the syringe with ready-to-use disinfectant wipes (e.g., CavWipes™ based on the active substances alcohols and quaternary compounds (QAV)) or equivalent disinfection agents.
- Follow the disinfection agent manufacturer’s directions for use, observing in particular the contact times.
- Disinfect the Penta Elastomer syringe for three minutes using at least one new disinfectant wipe. Be sure that all of the surfaces of the Penta Elastomer syringe remain moist during the entire three-minute period.

Manual Cleaning and Disinfection of the Knurled Nut, Plunger, and Handle of the Penta Elastomer Syringe

An ultrasonic bath or, alternatively, a cleaning bath may be used for cleaning. Disinfect the syringe in a disinfection bath after cleaning.

Cleaning in an Ultrasonic Bath

- Cleaning agent: Use for example Neodisher® MediClean Forte, Dr. Weigert, 0.5% concentration, or an equivalent cleaning agent for at least 10 minutes.
- Rinse off the gross impurities under running water (use potable water or drinking water quality as a minimum) for 1 minute.
- Follow manufacturer’s instructions to prepare a cleaning solution in the ultrasonic bath.
- Place the knurled nut, plunger, and handle of the Penta Elastomer syringe in a wire basket, immerse it completely in the cleaning solution, and treat with ultrasound for 10 minutes.

- When no more visible impurities can be seen, use potable water or water of drinking water quality as a minimum to rinse off the cleaning solution residue for 1 minute.
- Use running, clean cloth to dry the knurled nut, plunger, and handle of the Penta Elastomer syringe.

Manual Cleaning in a Cleaning Bath

► Cleaning agent: Use for example Neodisher® MediClean Forte, Dr. Weigert, 0.5% concentration, or an equivalent cleaning agent for a minimum exposure time of 10 minutes.

- Place the knurled nut, plunger, and handle of the Penta Elastomer syringe in the cleaning bath for the specified exposure time. Be sure that they are adequately covered in the solution and do not touch each other. Enhance the cleaning by completely brushing all of the inner and outer surfaces with a soft brush.
- Remove the parts from the cleaning bath and rinse thoroughly at least three times for a minimum of 1 minute with water (of drinking water quality as a minimum).
- Check the parts after cleaning (see “Inspection, Maintenance, Testing”).

Manual Disinfection (High Level Disinfection)

- Disinfection agent: Use for example Cidex OPA®, Johnson & Johnson, undiluted, (or an equivalent disinfection agent) for a exposure time of 5 minutes.
- Place the cleaned and inspected knurled nut, plunger, and handle of the Penta Elastomer syringe in the disinfection bath for the specified exposure time (> 20° C). Be sure that the syringe parts are adequately covered in the solution and do not touch each other.

- Remove the knurled nut, plunger, and handle of the Penta Elastomer syringe from the disinfection bath and rinse thoroughly with water (of drinking water quality as a minimum) at least five times for a minimum of 1 minute.

Drying

Dry the knurled nut, plunger, and handle of the Penta Elastomer syringe at a clean place.

Inspection, Maintenance, Testing

Store the Penta Elastomer syringe under dry and dust-free conditions.

Check the knurled nut, plunger, and handle of the Penta Elastomer syringe for damage, discoloration, and contamination before every use and discard damaged parts.

If there is still contamination on any of the parts, they must be cleaned and disinfected again.

Notes

- Please report a serious incident occurring in relation to the device to 3M and the local competent authority (EU) or local regulatory authority.
- Penta Elastomer Syringe Tips are single-use devices and would be damaged when removing set impression material.
- Please keep the single-use accessory in the original bag to guarantee the traceability of the batch.

Disposal

Dispose of the contents or container in accordance with the applicable regulations. Please pay special attention to the disposal of contaminated waste to avoid health risks because of improper handling.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M Deutschland GmbH warrants this product will be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user’s application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH’s sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability

Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence, or strict liability.

Symbol Glossary

| Reference Number and Symbol Title | Symbol | Description of Symbol |
|--|---|---|
| ISO 15223-1 5.1.1 <p>Manufacturer</p> |  | Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. |
| ISO 15223-1 5.1.3 <p>Date of Manufacture</p> |  | Indicates the date when the medical device was manufactured. |
| ISO 15223-1 5.1.4 <p>Use-by date</p> |  | Indicates the date after which the medical device is not to be used. |
| ISO 15223-1 5.1.5 <p>Batch code</p> |  | Indicates the manufacturer’s batch code so that the batch or lot can be identified. |
| ISO 15223-1 5.1.6 <p>Catalogue number</p> |  | Indicates the manufacturer’s catalogue number so that the medical device can be identified. |
| ISO 15223-1 5.4.2 <p>Do not re-use</p> |  | Indicates a medical device that is intended for one use or for use on a single patient during a single procedure. |
| ISO 15223-1 5.4.4 <p>Caution</p> |  | Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. |
| CE Mark |  | Indicates conformity to European Union Medical Device Regulation or Directive. |
| Medical Device |  | Indicates the item is a medical device. |
| Rx Only |  | Indicates that U.S. Federal Law restricts this device to sale by or on the order of a dental professional. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1). |
| Green Dot |  | Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe. |

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3M ESPE, Penta and Pentamix are trademarks of 3M or 3M Deutschland GmbH. Used under license in Canada.
3M, ESPE, Penta and Pentamix sind des marques de commerce de 3M ou 3M Deutschland GmbH. Utilisées sous licence au Canada. © 2019, 3M. All rights reserved.