# Duceram® Kiss NE-Bonder

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### Additions to the Instructions for Use for Duceram® Kiss

1. General instructions on processing the NE-Bonder for non-precious metals
For information on indications, contraindications, warning messages for medical devices, safety notes, interactions, technical specifications, transport and storage, liquid combination options, and explanations of symbols on the product labels, please consult the enclosed Quick Reference or Instructions for Use for Duceram® Kiss.

2. Indications
The NE-Bonder is a material that complements the Duceram® Kiss ceramic system. It is designed to optimize the adaptation of the ceramic veneer to non-precious metal frameworks.

3. Casting
When casting non-precious alloys, it is recommended to use only ceramic crucibles and to use only new material.

4. Framework preparation
During finishing, the framework must be freed of any sharp edges. The framework is then sandblasted using alumina 250 µm at a pressure of 3–4 bar. An oxide firing step is possible but not mandatory; it only serves to check the framework.

5. Application and processing
The NE-Bonder is mixed with Ducera® Liquid B or Ducera Liquid OCL universal and applied to the framework in a uniform semi-covering layer. Then the NE-Bonder is fired at 980°C, which ensures that it is sintered to the framework in a homogeneous fashion and the subsequent vitrification is performed in a uniform manner. This in turn ensures that the non-precious alloy bonds securely to the ceramic veneer.

| Firing | Pre-heating | Drying time | Heating rate | Final temperature | Holding time | Vacuum | Tempering | Stress-relief cooling |
|--------|-------------|-------------|--------------|-------------------|--------------|--------|-----------|-----------------------|
|        | °C          | min.        | °C/min.      | °C                | min.         | hPa    |           |                       |
| Bonder | 575         | 7:00        | 55           | 980               | 2:00         | 50     | –         | –                     |

The values listed here are intended for orientation only and should be regarded only as guidelines. Your firing results may differ. All firing results depend on the performance of the furnace used, which in turn depends on the make, model and age of the furnace. Therefore, the guideline values will have to be adapted individually for each firing. We recommend running a test firing cycle to evaluate the performance of the furnace used. We have compiled and checked all values and other data with great care. However, we cannot be liable for your results under any circumstances.

For further processing using the Duceram® Kiss veneering ceramic material, please follow the Instructions for Use and in particular the special firing parameters for non-precious alloys. To prevent discolouration, any exposed non-precious alloy connectors, cast abutments, etc.) should be polished with rubber tips or wheels or sandblasted after each firing step.

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