ENGLISH

General Information

3M[™] Filtek[™] Supreme Flowable Restorative is a low viscosity, visible-light activated, radiopaque, flowable nanocomposite. The restorable is packaged in syringes. It is available in a variety of tooth-colored shades. The shades offered with Filtek Supreme flowable restorative were designed to coordinate with shades offered with Filtek[™] Supreme Ultra/Supreme XTE/Ultimate/Z350 XT Universal Restorative, and Filtek[™] Universal Restorative. Filtek Supreme flowable restorative contains bisGMA, TEGDMA and Procyrlates resins. The fillers are a combination of ytterbium trifluoride filler with a range of particles sizes from 0.1 to 5.0 microns, a non-agglomerated/non-aggregated surface modified 20 nm silica filler, a non-agglomerated/non-aggregated surface modified 75 nm silica filler, and a surface modified aggregated zirconia/silica cluster filler (comprised of 20 nm silica and 4 to 11 nm zirconia particles). The aggregate has an average cluster particle size of 0.6 to 10 microns. The inorganic filler loading is approximately 65% by weight (46% by volume).

Indications

- Class III and V restorations
- Restoration of minimally invasive cavity preparations (including small, non-stress-bearing occlusal restorations)
- Base/liner under direct restorations
- Repair of small defects in esthetic indirect restorations
- Pit and fissure sealant
- Undercut blockout
- Repair of resin and acrylic temporary materials

Precautionary Information for Patients

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. Avoid use of this product in patients with known acrylate allergies. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

Precautionary Information for Dental Personnel

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

3M SDSs can be obtained from www.3M.com or contact your local subsidiary.

Instructions for Use

Preparation

1. **Prophy:** Teeth should be cleaned with pumice and water to remove surface stains.

2. **Shade Selection:** Prior to isolation of tooth, select the appropriate shade(s) of Filtek Supreme flowable restorative using a standard VITAPAN[®] classical shade guide. Shade selection accuracy may be enhanced via a mock-up; place the chosen shade of the restorative material on the unetched tooth. Manipulate the material to approximate the thickness and site of the restoration. Cure. Evaluate the shade match under different lighting sources. Remove the restorative material from the unetched tooth with an explorer. Repeat process until an acceptable shade match is achieved. For sealants, a contrasting shade may be desirable to enhance detection.

3. **Isolation:** A rubber dam is the preferred method of isolation. Cotton rolls and an evacuator can also be used.

Directions

1. **Pulp protection:** If a pulp exposure has occurred and the situation warrants a direct pulp capping procedure, use a minimum amount of calcium hydroxide on the exposure followed by an application of Vitrebond[™] or Vitrebond[™] Plus Light Cure Glass Ionomer Liner/Base, manufactured by 3M ESPE. Vitrebond or Vitrebond Plus liner/bases may also be used to line areas of deep cavity excavation.

2. **Adhesive System:** Follow manufacturer's instructions regarding placement.

Note: Silane treatment is recommended for repair of ceramic restorations followed by the adhesive application.

3. **Warning:** Syringe can be warmed up to 70°C for one hour up to 25 times in commercially available composite warmer devices.

Dispensing:

- 4.1. **Syringe delivery:** Filtek Supreme flowable restorative can be delivered directly from the dispensing tip.
 - 4.1.1. Protective eyewear for patients and staff is recommended when using the dispensing tip.
 - 4.1.2. Prepare the dispensing tip: Remove cap and save. After use, the dispensing tip need to be discarded and replaced with syringe storage cap.
 - 4.1.3. Twist the disposable dispensing tip securely onto the syringe. As needed, bend the metal cannula in any direction, up to a 90° angle to facilitate access to prep.
 - 4.1.4. To avoid contamination of the syringe during treatment, a standard dental protective sleeve should be used. See the section "Cleaning & Disinfection."
 - 4.1.5. Holding the tip away from the patient and any dental staff, express a small amount of Filtek Supreme flowable restorative to assure that the delivery system is not plugged.
 - 4.1.6. If plugged, remove the dispensing tip and express a small amount of Filtek Supreme flowable restorative directly from the syringe. Remove any visible plug, if present, from the syringe opening. Replace dispensing tip and again express composite. Filtek Supreme flowable restorative may be extruded onto a dispensing pad and applied with a brush or other appropriate instrument.

Placement:

- 5.1. Place and light cure restorative in increments as indicated in Section 6.
 - 5.2. Avoid intense light in the working field.
6. **Curing:** This product is intended to be cured by exposure to a halogen or LED light with a minimum intensity of 400 mW/cm² in the 400-500 nm range. Cure each increment by exposing its entire surface to a high intensity visible light source, such as a 3M ESPE curing light. Hold the light guide tip as close to the restorative as possible during light exposure.

Cure Time			
Shades	Increment Depth	All halogen lights (with output 400-1000 mW/cm ²)	Elipar [™] S10 and Elipar [™] Freelight 2 (LED lights with output 1000-2000 mW/cm ²)
Opaque Shades	1.5 mm	40 sec.	20 sec.
All Other Shades	2.0 mm	20 sec.	10 sec.

7. Complete the restoration:

- 7.1. **Base/Liner application:** Place a composite restorative material, such as Filtek Supreme Ultra/Supreme XTE/Ultimate/Z350 XT Universal Restorative directly over the cured Filtek Supreme flowable restorative. Follow the manufacturer's instructions regarding placement, curing, finishing, occlusal adjustment, and polishing.

7.2. Direct Restorative Application:

- 7.2.1. Contour restoration surfaces with fine finishing diamonds, burs or stones.
 - 7.2.2. Check occlusion with a thin articulating paper. Examine centric and lateral excursion contacts. Carefully adjust occlusion by removing material with a fine polishing diamond or stone.
 - 7.2.3. Polish with Sof-Lex[™] Finishing and Polishing System or with white stones, rubber points or polishing paste where discs are not suitable.
- 7.3. **Pit and Fissure Sealants:** Gently remove the inhibited layer remaining after light curing with slurry of pumice or polishing paste.

Storage and Use

1. This product is designed to be used at room temperature. If desired, the product may be warmed in a commercial composite warmer prior to use (at no higher than 70°C / 158°F, for no longer than 1 hour), up to 25 times.
2. If stored in cooler allow product to reach room temperature prior to use. Shelf life at room temperature is 24 months. Ambient temperatures routinely higher than 27°C/80°F may reduce shelf life. See outer package for expiration date.
3. Do not expose restorative materials to elevated temperatures, or to intense light.
4. Do not store materials in proximity to eugenol containing products.

5. Discard used syringe tip and replace with syringe storage cap.

Cleaning & Disinfection

1. Use new, unused disposable gloves during all clinical applications for the syringe to reduce the risk of cross-contamination.
2. Place the disposable tip on the syringe and place the protective sleeve over the syringe. Use the metal cannula to punch through the end of the protective sleeve.
3. Avoid contact between the reusable parts (e.g., the body of the syringe) and the patient's mouth.
4. After using the syringe in the protective sleeve, remove the sleeve carefully so that there is no contamination of the syringe from the outer surface of the protective sleeve.
5. After removing the protective sleeve and disposable tip, wipe the syringe with a ready-to-use cleaning cloth (e.g., CaviWipe[™]) for the recommended contact time on the wipe label.
6. Discard the used protective sleeve along with the syringe tip.
7. Always observe all applicable legal and hygiene regulations for dental offices and/or hospitals during the use and reprocessing of the device.

Note: Do not place the syringe into a disinfection bath or a cleaning-disinfection device (washer-disinfector).

Disposal – See the Safety Data Sheet (available at www.3M.com or through your local subsidiary) for disposal information.

Customer Information

No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty

3M warrants this product will be free from defects in material and manufacture. 3M MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and 3M's sole obligation shall be repair or replacement of the 3M product.

Limitation of Liability

Except where prohibited by law, 3M will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
CE Mark 0123		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
Rx Only		Indicates that the U.S. Federal Law restricts this device to sale by or on the order of [specific licensed provider]. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
CE Mark 0123		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
Rx Only		Indicates that the U.S. Federal Law restricts this device to sale by or on the order of [specific licensed provider]. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
CE Mark 0123		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
Rx Only		Indicates that the U.S. Federal Law restricts this device to sale by or on the order of [specific licensed provider]. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
CE Mark 0123		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
Rx Only		Indicates that the U.S. Federal Law restricts this device to sale by or on the order of [specific licensed provider]. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)

Symbol Glossary		
Symbol Title	Symbol	Description and Reference
Manufacturer		Indicates the medical device manufacturer as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC. Source: ISO 15223, 5.1.1
Authorized Representative in European Community		Indicates the authorized representative in the European Community. Source: ISO 15223, 5.1.2, 2014/35/EU, and/or 2014/30/EU
Date of Manufacture		Indicates the date when the medical device was manufactured. Source: ISO 15223, 5.1.3
Use-by date		Indicates the date after which the medical device is not to be used. Source: ISO 15223, 5.1.4
Batch code		Indicates the manufacturer's batch code so that the batch or lot can be identified. ISO 15223, 5.1.5
Catalogue number		Indicates the manufacturer's catalogue number so that the medical device can be identified. Source: ISO 15223, 5.1.6
Temperature limit		Indicates the temperature limits to which the medical device can be safely exposed. Source: ISO 15223, 5.3.7
Caution		Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot, for a variety of reasons, be presented on the medical device itself. Source: ISO 15223, 5.4.4
Unique device identifier		Indicates bar code to scan product information into patient electronic health record.
Importer		Indicates the entity importing the medical device into the EU.
CE Mark 0123		Indicates conformity to all applicable European Union Medical Device Regulations or and Directives with notified body involvement.
Rx Only		Indicates that the U.S. Federal Law restricts this device to sale by or on the order of [specific licensed provider]. 21 Code of Federal Regulations (CFR) sec. 801.109(b)(1)
Green Dot		Indicates a financial contribution to national packaging recovery company per European Directive No. 94/62 and corresponding national law. Packaging Recovery Organization Europe.
Low-density polyethylene		Indicates plastic component made of low-density polyethylene is recyclable. Official Journal of the EC; Commission Decision (97/129/EC)

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For more information see, HCBRegulatory.3M.com

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