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## GB Processing instructions

### 1 Intended purpose

Denture base materials are intended for the fabrication of splints, partial dentures, full dentures and model dentures.

### 2 Description of product and users

#### 2.1 Product description

Denta-Don is a heat-curing acrylic for the production of acrylic dentures using the injection or flask technique. The large range of colours allows the form and colour to be imitated perfectly.

#### 2.2 Users

For use by laboratory technicians in a dental laboratory

### 3 Composition

**Powder:** Pearl polymerisable made from polymethyl methacrylate (PMMA) pigments, initiators

**Liquid:** Methyl methacrylate, tetramethylene dimethacrylate, initiators, stabilisers

### 4 Indications

- For the production of partial and full dentures using the flask or injection technique
- For the production of splints
  - o Mouthguards
  - o Bite splints

### 5 Contraindications

- Do not use in the case of a known allergy to one of the components.

### 6 Warnings

**Liquid:** Contains methyl methacrylate, tetramethylene dimethacrylate

Danger. Highly flammable liquid and vapour. Causes skin irritation. May cause an allergic skin reaction. May cause respiratory irritation.

**Powder:** Contains dibenzoyl peroxide. May produce an allergic reaction. Harmful to aquatic life with long lasting effects.

### 7 Safety instructions

**Liquid:** Keep away from heat, hot surfaces, sparks, open flames and other ignition sources. No smoking. Use explosion-proof [electrical/ventilating/lighting] equipment. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water [or shower]. Store locked up.

**Powder:** Avoid release to the environment.

### 8 Interactions with other medicinal products

The dentist should consider known interactions and cross-reactions of the medical device with other materials already in the patient's mouth before using the product.

### 9 Application / Preparation

#### 9.1 Production of acrylic dentures

##### 9.1.1 Preparatory work

Produce the plaster cast and wax-up of the teeth as usual.

The fabricated tooth set in wax is secured with plaster or silicone (matrix/flask). Then the wax is eliminated as usual.

To prevent artificial teeth from breaking off, they must be roughened, have undercuts made, and be wetted with monomer or a suitable bonder.

The surfaces of the plaster models are separated

with an alginate separating agent according to the manufacturer's instructions.

#### 9.1.2 Mixing ratio

• 2.5 parts by weight powder and 1 part by weight liquid

• If a longer processing window is desired, we recommend a mixing ratio of 2.0 parts by weight powder and 1.0 part by weight liquid.

• Veins can be included, if desired. Please observe separate manufacturer's information and instructions.

#### 9.1.3 Processing in the pressing technique

After adding the powder to the liquid, mix Acryl BH thoroughly with a spatula in a mixing vessel for 30 seconds to form a homogeneous mixture.

The acrylic dough is workable after it has been allowed to thicken for approx. 10 min.

The processing window at approx. 21°C is:

- approx. 1 hour when mixed 2.5 : 1;
- approx. 2 hours when mixed 2.0 : 1.

After thickening, the acrylic dough must be kneaded for approx. 30 sec.

When the Denta-Don becomes viscoelastic, pack it into the prepared flask and press.

Maintain the pressure for approx. 15 minutes.

Secure the flask in a clamp for polymerisation.

Polymerise Denta-Don in the closed flask with a clamp for 20 min. in boiling water (or in a pressure vessel at 95°C). Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

#### 9.1.4 Processing in the UNIPRESS system (Brand: Schütz Dental. observe separate instructions):

Flasks and plaster surfaces are prepared as described in the system instructions. Denta-Don is mixed as described in the pressing technique and placed in the cylinder and injected. The polymerisation process corresponds to the pressing technique.

#### 9.1.5 Finishing the material

The dental appliance is finished using suitable instruments.

Sandpaper or pumice stone is used for pre-polishing.

### 9.2 Production of splints

#### 9.2.1 Preparatory work

An impression is taken and a plaster model (cast) produced. Undercuts are blocked out with suitable materials. The splint is then modelled in wax as desired. A plaster index is created and the wax removed.

The plaster model is separated using suitable materials according to manufacturer's instructions.

#### 9.2.2 Mixing ratio

• 2.5 parts by weight powder and 1 part by weight liquid

• If a longer processing window is desired, we recommend a mixing ratio of 2.0 parts by weight powder and 1.0 part by weight liquid.

#### 9.2.3 Processing

After adding the powder to the liquid, mix Acryl BH thoroughly with a spatula in a mixing vessel for 30 seconds to form a homogeneous mixture.

The acrylic dough is workable after it has been allowed to thicken for 10 min.

Model acrylic on the model, close flask. Secure in clamp and press.

Maintain the pressure for approx. 15 minutes.

Secure the flask in a clamp for polymerisation.

Polymerise Denta-Don in the closed flask with a clamp for 20 min. in boiling water (or in a pressure vessel at 95°C). Allow the flask to cool to room temperature in the clamp and then deflask in the usual way.

#### 9.2.4 Finishing the material

After polymerisation, the splint is lifted from the model and finished using suitable instruments.

Silicone polishers (square-edge or knife-edge wheels or cylinders), tungsten carbide burs and diamond instruments are suitable for finishing.

Sandpaper or pumice stone is used for pre-polishing.

#### 9.3 Polishing

The material is polished with goat bristle brushes, polishing paste and soft cotton wool discs. Careful surface finishing and polishing is necessary for an optimal result.

#### 9.4 Repairs

Repairs can be made with cold-curing polymers according to manufacturer's instructions. Before applying the acrylic dough, the fractured surfaces

must be roughened, cleaned and wetted with monomer. The polymerisation time in a pressure polymerisation device is 15 min., at a pressure of 2-5 bar and a water temperature of 45°C. Please observe separate manufacturer's instructions.

### 10 Troubleshooting / FAQ List\*

#### 11 Storage and handling information

Storage temperature 10-25°C / 50-77°F. Close container carefully.

#### 12 Shelf life

The maximum shelf life is printed on the label of each pack. Do not use after the expiry date.

#### 13 Warnings on side effects

With proper preparation and use of this medical device, adverse effects are extremely rare. However, immune reactions (such as allergies) or local discomfort cannot in principle be ruled out completely. All serious incidents which occur in connection with the use of this product are to be reported to the manufacturer indicated below and the competent authority in each case.

#### 14 Instructions for disposal

Leftover quantities and packaging materials are to be disposed of according to the local and/or statutory regulations.

### \* Troubleshooting

Fault	Cause	Corrective action
Acrylic does not cure properly	Incorrect mixing ratio	– Refer to the instructions for use – Ensure correct mixing ratio
Acrylic is whitish at base	– Inadequate insulation – Cast not steeped in water	– Ensure adequate insulation – Check insulating materials if necessary – Allow cast to steep for approx. 10 min in lukewarm water
Acrylic forms bubbles at base	– Cast not steeped in water	Allow cast to steep for approx. 10 min in lukewarm water
Acrylic surface white/porous	Acrylic injected too late or placed in pressure vessel too late	Observe the processing times
Differences in colouration following repairs	Incorrect mixing ratio	Ensure correct mixing ratio
Acrylic breaks when lifted out	Inadequate insulation	Check insulation
Acrylic/denture has whitish streaks	– Inadequate mixing – Too much powder used	– Carefully stir acrylic dough with spatula – Ensure correct mixing ratio
Splint breaks	Contouring too thin/unsuitable	Ensure sufficient layer thickness (at least 2-3 mm)

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