PERMABOND PINK / DENTIN ADHÄSION

Stand: 07/2018

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Stand:070618

GB Instructions for use

If there is anything in this instruction for use that you do not understand, please contact our customer service department before using this product.

Product description:

M+W Permabond pink/ dentin Adhäsion is a light-curing paste opaquer for clasp coating and the coating of model cast retentions with an integrated adhesive component.

Thus, it is not necessary to use an additional adhesive bonding system.

Application:

After the casting, the clasp / frame is processed as usual. The surface is sand-blasted with 50 µm aluminum oxide (2 bar) and blast-cleaned with oil-free compressed air. A thin layer of M+W Permabond pink/ dentin Adhäsion is applied with a stiff shorthair brush. To ensure that the coating covers the whole surface with the required thickness, at least two layers must be applied. A transparent opaquer layer has a negative effect on the color. The flowability of the opaquer can be optimized by stirring it on a mixing pad. In the area of retentions, make sure that the applied layer is extra thin. Each layer is polymerized separately.

Spektra LED:	30 sec.
Spektra 2000:	3 min.
PLC Spektra:	3 min.

Opaquer layers that have been applied too thickly prevent an optimal polymerization and thus weaken the bond. The hardening cannot be improved by extending the polymerization time.

Sources of error / Error prevention:

Remove any layers of grease and polishing waste from the frame surfaces. Only apply very thin layers of opaquer. Do not mix M+W Permabond pink/ dentin Adhäsion with other opaquer liquids / powders. Do not use intermediate polymerization units for polymerization.

Processing of the materials:

Depending on lighting conditions 1 - 3 min.

Storage:

The material should be stored at 10 - 25° C. Close the syringe immediately after use and with due diligence and protect it against direct sunlight.

Best before date:

The best before date is printed onto the label of each syringe. Do not use after this date.

Syringe content:

2 g – Permabond dentin Adhäsion

5g – Permabond pink

Guarantee

Our technical instructions, regardless of whether they are provided verbally, in writing or during practical demonstrations, are based on our own experience and should only be considered guidelines. As our products are subject to continued development, we reserve the right to modify them.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us.

Contra-indications / interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please supply the dentist with the above information, if this medical device is used to produce a special model. Please also consider the safety data sheets.

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany · Tel.: +49 (0)6003 814-0 · Fax: +49 (0)6003 814-901