GLASFILL

Gebräuchsanweisung - Iso 9917-1:2007 Glas-Polyalkenoat Klasse 4.2c

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heat and desiccation of the cement. After finishing, the surface should be coated with varnish for protection of moisture.

WARNING:

Do not expose patients or users known to be allergic to this type of material. Avoid contact of liquid and powder with oral mucosa, eyes, and skin. In case of such contact, wash thoroughly with large quantities of water and obtain medical advice as soon as possible.

STORAGE:

Store in a cool, dry place. The powder is protected against moisture uptake by a desiccant capsule, but the cap should be replaced immediately after use.

BATCH CODE:

[LOT] The batch code gives an open date of manufacture in month, year, day format with a numerical suffix to uniquely identify the batch of material. Please quote this batch number in all correspondence.

⌚

The expiry date is shown in year, month format. Do not use the product after this date. For test purposes, the powder-liquid ratio is 1.05 g to 0.15 g at 23±1°C and 50±10% RH. ISO 9917 net Setting Time:

Express: 2 to 3 minutes from end of mix at 36-38°C.

This product has been formulated for use in dentistry and is intended for use by dental professionals only.

Directions for use - ISO 9917-1:2007 Glass Polyalkenoate class 4.2c

Glass Ionomer restorative cement powder is a fast setting glass polyalkenoate material activated by mixing with water. The mixed cement adheres directly to enamel, dentine and cementum. Contains fluoride.

INDICATIONS:

- 1) Class V cavities (particularly erosion lesions).
- 2) Class III cavities.
- 3) Class I and II cavities in deciduous teeth.
- 4) Pits and fissures.

CONTRA-INDICATIONS:

Pulp capping.

CONTENTS

Powder 10g, dispensing bottle, dropper tip, measuring scoop, instructions for use.

PROCEDURE / PRECAUTION

(1) CAVITY PREPARATION:

Glass Ionomer restoratives, because of their adhesive qualities, can be used with minimal cavity preparation techniques. (Adhesion to set amalgam restorations may be useful in certain situations.) Cavity liner is required only in the deepest cavities where calcium hydroxide should be used. In class V erosion lesions, 25% polyacrylic acid solution should be applied on a pledget of cotton wool for 30 seconds to clean the tooth surface. Wash with water and dry using oil-free air. In areas which are heavily coated with plaque or tartar, prophylaxis using oil-free prophylaxis paste prior to tooth cleansing is required. In cut cavities, the use of polyacrylic acid solution is not required unless subsequent contamination with saliva occurs.

(2) MIXING:

Mix ratio 2 scoops of powder:2 drops of water (powder: liquid ratio 7:1 m:m at 22-24°C and 40-60% RH). Invert bottle to fluff powder; this ensures correct scoop weight. Do not compress powder against side of the bottle with the scoop. Remove excess powder using a flat spatula blade, again taking care to avoid compressing the powder. To deliver accurate water drops, hold bottle vertically. If water does not drop from the TIP of the steel tube, clean this with a tissue. Use a glass block for best results and a stainless steel spatula. Incorporate half the powder into the water as quickly as possible (5-10 seconds) and then add the remainder and spatulate to a THICK putty-like consistency. Total mixing time 30 seconds. DO NOT ADD POWDER IN SMALL INCREMENTS.

(3) PLACEMENT:

Apply to cavity using normal instruments (non-stick aluminium instruments are particularly suitable). To avoid cement adhering to steel instruments, dip the clean instrument either into the powder or methylated spirit. If stainless steel matrices are used, they should be lightly coated with petroleum jelly.

(4) WORKING TIME:

Express Set: Minimum 1 minute 20 seconds from start of mix at 22-24°C.

(5) CLINICAL SETTING TIME:

Express Set: Maximum 3 minutes from end of mix.

(6) FINISHING:

Best results are obtained by polishing with abrasive discs and stones using water spray lubrication during a subsequent visit. At 5 minutes after placement the material is hard enough to finish using abrasive discs and stones, but petroleum jelly must be used as lubricant to prevent excess

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