ENGLISH
<b>General Information</b>
Vitrebond™ Plus light cure glass ionomer manufactured by 3M ESPE is a two-part liquid/paste system. The liquid/paste materials are contained in the Clicker™ Dispensing System manufactured by 3M ESPE. This dispensing system provides simultaneous dispensing of each component for a consistent mix.
The composition is based on Vitrebond™ light-cure glass ionomer liner/ base. The paste contains a radiopaque fluoro-aluminosilicate glass. The liquid contains a modified polyalkylene acid. Vitrebond Plus liner/base provides the major benefits of glass ionomer cements including adhesion to wet structure and sustained fluoride release. Additionally, Vitrebond Plus liner/base offers a combination of a prolonged working time with a short set time achieved by light curing.

Indications
Vitrebond Plus liner/base is indicated for lining and basing applications under the following restorations:

- Composite
- Amalgam
- Ceramic
- Metal

#### Contraindications:

Vitrebond Plus liner/base is **not indicated** for **direct pulp capping**. If a pulp exposure has occurred and if the situation warrants a direct pulp capping procedure use a minimum amount of calcium hydroxide on the exposure followed by an application of Vitrebond Plus liner/base.

#### Precautionary Information for Patients:

Avoid use of this product in patients with known acrylate allergies. This product contains substances that may cause an allergic reaction by skin contact in certain individuals. If prolonged contact with oral soft tissue occurs, flush with large amounts of water. If allergic reaction occurs, seek medical attention as needed, remove the product if necessary and discontinue future use of the product.

#### Precautionary Information for Dental Personnel:

This product contains substances that may cause an allergic reaction by skin contact in certain individuals. To reduce the risk of allergic response, minimize exposure to these materials. In particular, avoid exposure to uncured product. If skin contact occurs, wash skin with soap and water. Use of protective gloves and a no-touch technique is recommended. Acrylates may penetrate commonly used gloves. If product contacts glove, remove and discard glove, wash hands immediately with soap and water and then re-glove. If allergic reaction occurs, seek medical attention as needed.

#### Instructions for Use

**1. Isolation:** A rubber dam is the preferred method of isolation. Avoid water and saliva contamination during application and setting of the Vitrebond Plus liner/base.

**2. Restoration/ooth Preparation:** Remove carious dentin and all amalgam or other base material from the internal form of the preparation. Rinse and dry cavity. Leave tooth surface moist.

**Do not overdry.**

**3. Pulp Protection:** Vitrebond Plus liner/base is **not indicated** for **direct pulp capping**. If a pulp exposure has occurred and if the situation warrants a direct pulp capping procedure use a minimum amount of calcium hydroxide on the exposure followed by an application of Vitrebond Plus liner/base.

**4. Dentin Pretreatment: Dentin pretreatment is not recommended.** The use of smear layer etchants such as polyacrylic acid based systems results in decreased adhesion of Vitrebond Plus liner/base.

#### 5. Dispensing

a) **Remove Cap:** Press and hold tab to unlock the protective cap. Slide cap off clicker dispenser.

b) **Dispense:** Touch dispenser tip to mix pad. Fully depress clicker lever to dispense "1 click" of Vitrebond Plus liner/base to the mix pad. Release lever when paste stops extruding (1-2 seconds). Repeat dispensing process for additional material. 1-2 clicks will be sufficient for most restorations.

c) **Clean:** Wipe dispenser tip with alcohol-dampened gauze.

d) **Replace Protective Cap:** Replace protective cap immediately after dispensing. Hold the sides of the clicker cartridge and side cap into place until securely latched as indicated by the audible "click".

**Do not depress the clicker lever during cap removal and/or replacement.**

**Do not advance the dispenser plunger during cap removal and/or replacement.**

An approximate number of Vitrebond Plus liner/base applications remaining in the Clicker Dispensing System are indicated on the underside of the plunger adjacent to the black arrow on the cartridge.

**6. Mixing:** Using a small spatula, mix paste/liquid components together rapidly (10-15 seconds). The mixed Vitrebond Plus liner/base should have a smooth consistency and glossy appearance.

To minimize water evaporation and maximize working time, confine spatulation of the paste/liquid to a small area on the mixing pad, about one inch (2.5 cm) in diameter.

**7. Application and curing:** Avoid water and saliva contamination during application and setting of the liner/base. Rubber dam is the best means of isolation.

Apply a thin layer (1/2mm or less) of the mixed Vitrebond Plus liner/ base material to the dentin surfaces of the prepared cavity using a ball applicator or other suitable instrument. Do not take out to the margins.

Vitrebond Plus liner/base material has a minimum working time of 2 minutes 30 seconds at a room temperature of about 73°F (23°C). Higher temperatures will shorten working time.

Cure Vitrebond Plus liner/base by exposing layers of 1.5 mm or less for 20 seconds to a 3M ESPE curing light manufactured by 3M ESPE or other curing unit of comparable intensity.

The delayed auto-setting mechanism of Vitrebond Plus liner/base will ensure an eventual cure. If material shielded from light polymerization such as undercut areas.

Where a thicker application of the Vitrebond Plus liner/base is desired, the best adhesion can be achieved by first placing and light curing a thin layer followed by placement of a second layer up to 1.5mm in depth and light cure for 20 seconds.

**8. Adhesive system:** Continue with the bonding step of the restorative procedure starting with etching. Etchant on the liner/base is not deleterious.

#### Storage and Use

1. The liner/base is designed for use at room temperature of approximately 10-27°C or 50-80°F. Avoid elevated temperatures.

2. Vitrebond Plus liner/base is light sensitive. Protect from ambient light exposure by dispensing just prior to use and replacing Clicker dispenser cap immediately after dispensing.

3. Clicker disinfection. Disinfect the capped Clicker using an intermediate level disinfection process (liquid contact) as recommended by the CDC (Centers for Disease Control) and endorsed by the ADA (American Dental Association). *Guidelines for Infection Control in Dental Health-Care Settings –2003* (Vol. 52; No. RR-17), Center for Disease Control and Prevention.

4. Shelf life of the Vitrebond Plus liner/base at room temperature is 24 months. See outer package for expiration dates.

5. Storage in refrigerator ensures longest possible shelf life. Allow to reach room temperature prior to use.

6. Do not store materials in proximity to eugenol containing products. No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

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