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PERMAPLAST BULK

GB INSTRUCTIONS FOR USE

If there is anything in this instruction for use that you do not understand, please contact our customer service department before using the product.

As the manufacturer of this medical device, we inform our users and patients that all serious events occurring in connection with it must be reported to us (the manufacturers) as well as the relevant authorities in the Member State where the user and/or patient is resident.

M+W Permaplast Bulk is a light curing, flowable, highly radio-opaque (210 % Al), composite for posterior restorations, for direct filling and the „bulk fill“ technique. M+W Permaplast Bulk is available in syringes and cavifils. The cavifils are for single use. Please do not reuse them, as this makes it impossible to rule out contamination and germ formation. The product complies with the guideline DIN EN ISO 4049.

Composition

Glass powder, diurethane dimethacrylate, tetramethylene dimethacrylate, silicon dioxide

Totale filler

77 % by weight (57 % by volume) inorganic filler (0,005-40µm)

Indications for use

- Fillings with layer thicknesses to 4 mm in Black's class I cavities.
- Base fillings/liners with layer thicknesses up to 4 mm in Black's class II cavities. For the following layers in a class II cavity, only use composites which are suited to filling Black's class II cavities.
- Cavity lining – as a first (bottom) layer in Black's class I and II cavities.

Application

Preparatory measures

Before beginning the treatment, clean the hard tooth substance with non-fluoride polishing paste.

1. Cavity preparation

Prepare the cavity with a minimally invasively method as generally required for adhesive techniques. All enamel margins in the anterior region must be bevelled. Do not bevel the margins in the posterior region and avoid slice preparations. Spray the cavity with water to clean it, remove all debris and dry it. Isolate the cavity to keep it dry. It is advisable to place a rubberdam.

2. Pulp protection / liner

The use of a liner can be foregone, if an enamel dentin adhesive is used. In case of very deep cavities that come close to the pulp, cover the

according surfaces with a calcium hydroxide compound.

3. Approximal contact areas

For cavities with approximal sections, place a transparent matrix and fix it in place.

4. Adhesive system

Etch and bond according to manufacturer's instructions.

5a. Application of M+W Permaplast Bulk Syringes

Lock the application tips on the syringe by turning it clockwise. Place layers of max. 4mm of M+W Permaplast Bulk directly into the cavity using the curved application tips supplied with the material. For hygienic reasons, these application tips must only be used once and then discarded. Avoid air bubbles when applying the material. Wet the prepared tooth surfaces thoroughly.

Application of M+W Permaplast Bulk Cavifils

Insert the cavifil into the dispenser. Remove the sealing cap. Place the cavifil in the correct angle towards the cavity. Inject the material into the cavity. Apply slow and steady pressure to the cavifil. Do not use excessive force! (max. layer thickness 4mm) To remove the cavifil from the dispenser after use, retract the plunger. Next, remove the cavifil. Please note: For hygiene reasons, cavifils are intended for single use only.

6. Polymerization

Polymerize each layer with a conventional polymerization unit for 20 seconds, holding the light guide as close as possible to the surface. During polymerization, a dispersion layer forms on the material surface. This layer must not be touched or removed. It forms the bond with the layer of composite applied subsequently (e. g. with a universal or posterior composite) or for creating aesthetic cusps. Alternatively, M+W Permaplast Bulk can be applied in a thin layer to form a liner on exposed dentin. The material can be polymerized in layers of up to 4 mm thickness.

7. Finishing

M+W Permaplast Bulk can be finished and polished directly after polymerization. To finish, use finishing diamond tools, flexible discs, silicone polishers and polishing brushes. Check occlusion and articulation. Grind down if necessary, thus avoiding premature or deflective contacts on the surface of the filling.

Please note

- When creating large and time-consuming restorations either cover the material with

a lightproof foil or temporarily move the operating light away from the site in order to avoid premature polymerization of the material.

- For hygienic reasons, the curved application tips supplied with the material must only be used once, as a contamination of the material and a transfer of germs cannot be ruled out.
- Use a polymerization unit with an emission spectrum of 350–500 nm to cure this material. As the required physical properties can only be achieved if the lamp is functioning correctly, its luminous intensity must be checked regularly as described by the manufacturer.

Light intensity for curing	>1200 mW/cm ²
Wavelength for curing	350-500 nm
Curing time	20 sec.

Hazard statements/precautionary statements

Contains Tetramethylene dimethacrylate

Warning: May cause an allergic skin reaction. Wear protective gloves/protective clothing/eye protection/face protection. If skin irritation or rash occurs: Get medical advice/attention.

Contra-indications

If a patient has known hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Side-effects

With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. Should you learn about unwanted side-effects – even if it is doubtful that the side-effect has been caused by our product – please kindly contact us. To prevent possible reactions of the pulp in cavities where the dentin is exposed, the pulp must be protected adequately (apply e. g. a calcium hydroxide preparation).

Interaction with other substances

As phenolic substances (such as eugenol) inhibit polymerization, do not use cavity liners containing such substances (e. g. zinc Eugenolcements).

Disinfection / Protection from cross-contamination

Place the functioning syringe with attached application tip into a suitably shaped barrier sheath; pierce end of sheath with cannula, exposing the cannula for use. Using a barrier sheath facilitates cleaning and disinfection

of the syringe between patients. After use of sheathed syringe, remove delivery tip and sheath by grasping on the hub of the delivery tip through the sheath, twist and remove tip along with sheath. Discard used tip and sheath in appropriate waste stream. Replace syringe storage cap. Disinfect – After removal of the delivery tip and sheath, disinfect this product using an intermediate level disinfection process (liquid contact) as recommended by the Centers for Disease Control and endorsed by the American Dental Association, Guidelines for Infection Control in Dental Health-Care Settings – 2003 (Vol. 52; No. RR-17), Centers for Disease Control and Prevention (USA).

Storage

Store at 10 – 25 °C (50 – 77 °F). Close syringes tightly immediately after use. Let the material reach room temperature before use. After use, retract the plunger of the syringe slightly to prevent the apertures from becoming blocked. Do not use after the expiry date (refer to label on syringe). For use by dentists only. Keep out of children's reach. This product was developed specifically for the described range of application. It must only be used according to the above instructions. The manufacturer will not be held liable for damages caused by incorrect handling or incorrect processing of the material.

Troubleshooting

Trouble	Cause	Remedy
Composite does not cure	Luminous intensity of the polymerization unit insufficient	Check luminous intensity; replace light source, if necessary
	Emitted spectral range of the polymerization unit insufficient	Consult manufacturer of polymerization unit; recommended spectral range: 350 – 500 nm
Composite seems to be too hard/firm inside the syringe/ cavifils	Material was stored at temperatures below 10°C for a longer period of time	Let composite reach room temperature before use
	Syringe was not closed tightly which caused part of the material to cure	close syringe tightly each time after removing composite
Composite does not cure sufficiently	Layer thickness per polymerization cycle too high	Keep to max. layer thickness of 4mm
Restoration seems too yellow when compared color reference	Insufficient polymerization of the composite layers	Repeat polymerization cycle several times, for a minimum of 20 sec

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GDF GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany · Tel.: +49 (0)6003 814-0 · Fax: +49 (0)6003 814-901
Vertrieb durch Müller & Weygandt GmbH · Reichardsweide 40 · D-63654 Büdingen · Tel.: +49 (0) 60 42 – 88 00 88