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## Intended use

Light-curing lacquer for plastic dental applications

## Indication for use

Light-curing lacquer for splints, partial and full dentures and temporary dental crowns in C&B technology.

## Indication

For smooth coating of plastic applications, to prevent mechanical irritation of the tongue and mucosa.

## Contraindication

Do not use in case of allergies to one of the ingredients.

## Patient target group

Patients with loss of teeth and dental substance defects, as well as well as patients who require splint therapy with plastic structures.

## Intended users

Dentists, dental technicians

## Requirements / equipment

|              |                                 |
|--------------|---------------------------------|
| EyeVolution® | FotoDent® flash (with nitrogen) |
|--------------|---------------------------------|

## Basic Material Properties

See technical datasheet (TDS). Technical datasheet available on request.

## Specific Manufacturing Considerations

### Caution

During processing, we recommend wearing personal protective equipment, e. g. protective work clothing, suitable gloves and safety goggles. Use a suitable extraction system or face mask during finishing.

## Process description

### Preparation

The area to be lacquered should be completely finished, dry and free from dust and grease residues. Due to the thin coat thickness, NanoVarnish is not suited as polish replacement. Therefore, the areas to be coated should be appropriately pre-polished.

## Lot number / Expiry date

The batch number and the use-by date are located on both the outer packaging and the product. When submitting a complaint about a product, please always quote the

## Further processing

For the best possible coating result, 10 MicroBrush are included in the package. These application aids ensure optimum distribution of even small amounts of material on the plastic surface.

The lacquer should only be applied quickly and as a thin, uniform film. Perform intermediate polymerization, if required. Only take the required quantity out of the bottle, and immediately close the bottle.

## Curing

|                                 |                       |
|---------------------------------|-----------------------|
| EyeVolution®                    | 7 min                 |
| FotoDent® flash (with nitrogen) | 1 min pre-gas + 4 min |

## Disposal

Disposal of contents / container in accordance with local / regional / national / international regulations and according to the safety data sheet.

## Hazards identification

Hazard and safety notes can be found on the product label and the corresponding safety data sheet.

## Cleaning

The plastic applications should be cleaned before first use in the patient using specifically available cleaning agents. Use of aggressive cleaning agents may damage the material and the surface.

## Disinfection

We recommend the use of a 0.2 % chlorhexidine solution or 0.5 % sodium hypochlorite solution with an exposure time of 5 minutes for disinfection, and then allow to dry completely.

## Sterilization

NanoVarnish cannot be sterilised.

## Serious incidents

All serious incidents occurring in relation to the product shall be reported to the manufacturer and to the competent authority of the Member State in which the user and/or patient is settled.

## Lot number / Expiry date

batch number. Do not use the product after expiration of the use-by date.

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