D

Warnhinweise

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GB

Instructions for preparation and sterilisation

Warnings

Instrumentes are not sterile upon receipt and must be sterilized before use in accordance with the following instruction. Instrumentes can be damaged by alkaline and acidic detergents. Always use a pH (7) neutral detergent or solution, Temperature do not exceed 134°C and remove all packaging prior to cleaning and sterilization, contaminated instrumentes should be handled with gloves and eye protection should be used.

Limitations

Repeated processing has minimal effect on these instrumentes. End of life is normally determined by wear and damage due to use. Never clean color-etched probes or black composite instrumentes in ultrasonic baths.

Preparation for Cleaning

No specific requirement, Disassemble not required. Rinsing: Immediately after surgery rinse instrumentes under warm (not hot) water. It may be helpful to use a nylon toothbrush to rinse the lock boxes and joints of the instrumentes. Be sure to remove all blood, body fluids, and tissue. Cleaning: If you do not clean your instrument immediately after rinsing, instrumentes should be submerged in a solution of water and neutral pH (7) detergent. They should never be placed in saline solution, as it may cause corrosion and eventually irreversible damage to the instrument.

Ultrasonic Cleaning

For micro and delicate instrumentes, use manual cleaning. Instrumentes should be processed in a cleaner for the full recommended cycle time—usually 5 to 10 minutes. Place instrumentes in open position into the ultrasonic cleaner. Make sure that “Sharp” instrumentes do not touch other instrumentes. All instrumentes have to be fully submerged. Do not place dissimilar metals (stainless, copper, chrome, plated, etc.) in the same cleaning cycle. Change solution freq.

Automated reprocessing

Simplee instrumentes are suitable for machine reprocessing in the thermal disinfector (WD). Please follow manufacturer’s instructions for the respective thermal disinfector (WD).

After Cleaning

If instrumentes are to be stored, let them, air dry and store them in a clean and dry environment.

Thermal disinfection

Instead of manual cleaning and disinfection, the

instrumentes can be reprocessed in an automated washer-disinfector (WD). Please follow the specific manufacturer’s instructions for the washer-disinfector.

Autoclaving

If instrumentes are to be reused or autoclaved maintain lubricate of all instrumentes which have any “metal to metal” action such as scissors, Hemostats, needle holder, self-retaining retractors, etc. Recommend surgical lubricants such as instrument milk are best. Do not use WD-40, oil or other industrial lubricants. Put instrumentes up for autoclaving either individually or in sets.

Inspection and Function Testing

Visually inspect to ensure all contamination has been removed. Check for distortion, damage and wear. Cutting edges should be free of defects. Discard damaged, worn or corroded instrumentes. Ensure that detachable tips are secure.

Drying

Dry using paper toweling or dry heat not exceeding 134 °C.

Packaging for Sterilization

Instrumentes may be loaded into dedicated instrumentes trays or general purpose sterilization trays, ensure that cutting edges are protected, For individual instrumentes disposable paper or plastic pouches are ideal. Make sure you use a wide enough pouch (at least 10 cm) for instrumentes with ratchet locks such as needle holders and hemostats so the instrumentes can be sterilized in an open (unlocked) position.

Sterilization

The instrumentes can be sterilized in designated container systems or general purpose sterilization containers individually packed in the standard polyethylene or Sterilization must be done by steam with fractionated pre-vacuum. Based on DIN EN 285 and DIN EN 554 is recommended at a temperature of 134 °C, a contact time of at least 5 minutes. Effective steam sterilization can be achieved with the following cycle: Temperature cycle type minimum exposure time Pre-vacuum to 134 °C - 137 °C for 5 minutes. To avoid damage the instrument, the temperature of 180 °C must not be exceeded.

Storage

Sterile, packaged instrumentes should be stored so that they are protected from dust, moisture, insects and vermin, and temperature extremes.

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