3M ESPE

RelyX™ Temp NE

- Ⓔ **Temporary Cement – Zinc Oxide Non-Eugenol**
- Ⓕ **Temporaire Zinkoxidzement – Eugenolfrei**
- Ⓖ **Ciment temporaire à base d'oxyde de zinc sans eugénol**
- Ⓗ **Cemento Provisorio a base di ossido di zinco senza eugenolo**
- Ⓘ **Cemento temporal de óxido de Zinco sin Eugenol**
- Ⓜ **Cimento Provisório de Óxido de Zinco sem Eugenol**
- Ⓝ **Tijdelijk non-eugenol zinkoxide cement**
- Ⓟ **Προσωρινή κονία Οξειδίου Ψευδούργου Ανευ Ευγενόλης**
- Ⓡ **Temporaire zinkoxidzement, eugenulfrei**
- Ⓢ **Väläikaisemistis Sinkioksiidi, eugenofiton**
- Ⓣ **Provisorisk zinkoxid cement, eugenofrit**
- Ⓥ **Midertidig zinkoxidbaseret sement uten eugenol**

Instructions for Use
Gebruiksaanwijzing
Информация по Употреблению
Informação per Uso
Información de empleo
Informação de uso

Gebruiksaanwijzing
Οδηγίες χρήσης
Brúksanvisning
Κατάχρηση
Bugsanvisning
Brúksanvisning

3M Deutschland GmbH
 Dental Products
 Carl-Schubert-Str. 1
 41460 Neuss - Germany

3M ESPE Dental
 3M, Paul. MN. 55144-1000 U.S.A.

3M ESPE Customer Care/USIS Information
 U.S. 1-800-426-2298 and Canada 1-888-363-3665.

3M ESPE and RelyX are trademarks of 3M Company.
 3M Deutschland GmbH, Ltd. under license in Canada.
 © 2012, 3M. All rights reserved.

ENGLISH
Description of Product
 RelyX™ Temp NE is an eugenol-free, self-curing cement to be used for temporary fixation of dental restorations. By adding Vaseline, the final hardness of the cement can be used for easier removal. In particular, RelyX Temp NE does not affect the polymerization of filling and luting composites. RelyX Temp NE is especially suitable for use on patients with sensitivity to eugenol.

* These Instructions for Use must be kept for the duration of product use.

Fields of Application
 Temporary fixation of provisional restorations,
 Temporary fixation of crowns, bridges, inlays, and onlays,
 Fixation of crowns and bridges on implant abutments

Precautionary Measures
For Patients
 The catalyst paste contains substances that may cause an allergic reaction in case of skin contact in certain individuals. Avoid use of this product in patients with known allergies to cadherins and its derivatives.
 If prolonged contact with oral tissue occurs, flush with large amounts of water. In an allergic reaction occurs, medical attention is needed, remove the product as soon as possible.

For Dental Personnel
 The catalyst paste contains substances that may cause an allergic reaction in case of skin contact in certain individuals. To reduce the risk of an allergic response, avoid exposure to this material. In particular, do not inhale the product. If skin contact occurs, wash skin with soap and water.
 The use of protective gloves and a no-touch technique is recommended. If the hands are soiled, remove gloves, and then re-glove. If an allergic reaction occurs, seek medical attention as needed.
 3M MSDS can be obtained from www.3m.com, or contact your local subsidiary.

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Application
 Before cementing, clean the cavity surfaces with a water spray and then dry. Apply a thin layer of resin on the surface to be cemented, and place the catalyst pastes with slight pressure.

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	min:sec	with maximum amount of Vaseline min:sec
Mixing	0:30	0:30
Working after start of mixing	1:30	2:30
Setting after start of main	3:30	4:00

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3M MSDS can be obtained from www.3m.com, or contact your local subsidiary.

Prevention of the Pulp
 In case of deep cavities, apply the protection of the pulp usually employed for temporary fixation.

Dispensing and Mixing
 Use RelyX Temp NE at an ambient temperature of 18° to 23°C/64° to 73°F and relative humidity of 45% to 55%.
 Dispense two equal lengths of base paste and catalyst onto a mixing pad.
 Dispensing the catalyst will increase the bonding strength and prolong the working time.
 Close the tubes immediately after use. Avoid interchanging the tube caps to prevent the possible hardening of the paste in the tube as a result of an unintentional carry-over of the pastes.
 The final hardness of the cement can be reduced by adding Vaseline in order to increase easier removal, e.g., for final bridges. This simultaneously extends the working time.
 Do not use a maximum of 42 mg of paste (a total of Vaseline if a length of 1 cm of base paste (172 mg) and a catalyst (78 mg) is used).
 Do not use plastic spatulas to thoroughly mix the pastes for 30 seconds until a homogeneous consistency is achieved.

Application
 Before cementing, clean the cavity surfaces with a water spray and then dry. Apply a thin layer of resin on the surface to be cemented, and place the catalyst pastes with light pressure.

Times
 The following times apply for an ambient temperature of 23°C/73°F:

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Setting is accelerated over an ambient temperature of 23°C/73°F and relative humidity over 55%.

Removal of Excessive Material
 Use an instrument to remove any excessive material after setting has occurred.

After Cementation
 Thoroughly cleanse the surface of the treated teeth and surrounding dentition with any residual cement from the matrix.

Storage and Shelf-Life
 Store the product at 15°–25°/59°–77°F.
 Do not use after the expiry date.

Customer Information
 No person is authorized to provide any information which deviates from the information provided in this instruction sheet.

Warranty
 3M Deutschland GmbH warrants this product to be free from defects in material and manufacture. 3M Deutschland GmbH MAKES NO OTHER WARRANTIES INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. User is responsible for determining the suitability of the product for your application. If this product is defective within the warranty period, your exclusive remedy and 3M Deutschland GmbH's sole obligation shall be repair or replacement of the 3M Deutschland GmbH product.

Limitation of Liability
 Except where prohibited by law, 3M Deutschland GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.
 Information as of: December 2011

Information for Patients
 The catalyst paste contains substances that may cause an allergic reaction in case of skin contact in certain individuals. To reduce the risk of an allergic response, avoid exposure to this material. In particular, do not inhale the product. If skin contact occurs, wash skin with soap and water.
 The use of protective gloves and a no-touch technique is recommended. If the hands are soiled, remove gloves, and then re-glove. If an allergic reaction occurs, seek medical attention as needed.
 3M MSDS can be obtained from www.3m.com, or contact your local subsidiary.

Information for Dental Personnel
 The catalyst paste contains substances that may cause an allergic reaction in case of skin contact in certain individuals. To reduce the risk of an allergic response, avoid exposure to this material. In particular, do not inhale the product. If skin contact occurs, wash skin with soap and water.
 The use of protective gloves and a no-touch technique is recommended. If the hands are soiled, remove gloves, and then re-glove. If an allergic reaction occurs, seek medical attention as needed.
 3M MSDS can be obtained from www.3m.com, or contact your local subsidiary.

Prevention of the Pulp
 In case of deep cavities, apply the protection of the pulp usually employed for temporary fixation.

Dispensing and Mixing
 Use RelyX Temp NE at an ambient temperature of 18°–23°C/64°–73°F and relative humidity of 45% to 55%.
 Dispense two equal lengths of base paste and catalyst onto a mixing pad.
 Dispensing the catalyst will increase the bonding strength and prolong the working time.
 Close the tubes immediately after use. Avoid interchanging the tube caps to prevent the possible hardening of the paste in the tube as a result of an unintentional carry-over of the pastes.
 The final hardness of the cement can be reduced by adding Vaseline in order to increase easier removal, e.g., for final bridges. This simultaneously extends the working time.
 Do not use a maximum of 42 mg of paste (a total of Vaseline if a length of 1 cm of base paste (172 mg) and a catalyst (78 mg) is used).
 Do not use plastic spatulas to thoroughly mix the pastes for 30 seconds until a homogeneous consistency is achieved.

Application
 Before cementing, clean the cavity surfaces with a water spray and then dry. Apply a thin layer of resin on the surface to be cemented, and place the catalyst pastes with light pressure.

Times
 The following times apply for an ambient temperature of 23°C/73°F:

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Dispensing and Mixing
 Use RelyX Temp NE at an ambient temperature of 18° to 23°C/64°–73°F and relative humidity of 45% to 55%.
 Dispense two equal lengths of base paste and catalyst onto a mixing pad.
 Dispensing the catalyst will increase the bonding strength and prolong the working time.
 Close the tubes immediately after use. Avoid interchanging the tube caps to prevent the possible hardening of the paste in the tube as a result of an unintentional carry-over of the pastes.
 The final hardness of the cement can be reduced by adding Vaseline in order to increase easier removal, e.g., for final bridges. This simultaneously extends the working time.
 Do not use a maximum of 42 mg of paste (a total of Vaseline if a length of 1 cm of base paste (172 mg) and a catalyst (78 mg) is used).
 Do not use plastic spatulas to thoroughly mix the pastes for 30 seconds until a homogeneous consistency is achieved.

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