Nahtmaterial Nylon

Stand: 170624

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Ⓤ Instructions for use

This package insert is not a reference to surgical techniques. It is designed to assist in using this product.

Product description

Nylon is a monofilament non-absorbable surgical suture material made of polyamide 6.0 or 6.6. Owing to its uniformly smooth surface nylon slips easily through tissue. Nylon is particularly noted for its very good knotting features and high tensile strength. Nylon is dyed in blue or black.

Nylon is available in a range of gauge sizes and lengths, nonneedled or attached to atraumatic stainless steel needles of varying types. The sutures are individually packed in sterile foil packs.

Nylon complies with the requirements of the European Pharmacopoeia and the United States Pharmacopoeia.

Indications

Nylon is a monofilament non-absorbable surgical suture material made of polyamide 6.0 or 6.6. Owing to its uniformly smooth surface nylon slips easily through tissue. Nylon is particularly noted for its very good knotting features and high tensile strength. Nylon is dyed in blue or black.

Nylon is available in a range of gauge sizes and lengths, nonneedled or attached to atraumatic stainless steel needles of varying types. The sutures are individually packed in sterile foil packs.

Nylon complies with the requirements of the European Pharmacopoeia and the United States Pharmacopoeia.

Action

This suture materials trigger a slight acute inflammatory reaction in the tissues. This leads to a general encapsulation of the suture material by fibrous connective tissue. Surgical suture material made of nylon is non-absorbable. Advanced hydrolysis of the nylon thread in vivo can lead to a decrease in the tensile strength of the thread.

Contra-indications

Owing to the in vivo decreasing tensile strength of surgical suture material made of nylon, surgical suture material made of nylon should not be used in patients, if a permanent retention of the tensile strength is necessary.

These sutures should not be used:

- in cardiovascular and neurological tissues and microsurgery.
- The use of this suture can be improper in elders, patients that are not properly fed, weak patients, and patients under certain conditions that could delay the healing of the wound.

Adverse reactions

Adverse reactions associated with this device include: wound dehiscence, failure to provide adequate wound support in sites where expansion, stretching or distension occur or in patients with delayed wound healing, transient local irritation at the wound site, transient inflammatory foreign body response, induration, erythema at the wound side and calculi formation. Like all foreign bodies these sutures may potentiate an existing infection.

Handling Instructions / Warnings / Precautions

Users should be familiar with surgical procedures and techniques and tie techniques involving non-absorbable sutures before employing these sutures for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used.

Particular factors of each patient should be considered regarding healing process in vivo.

Skin, vaginal and conjunctival sutures which must remain in place longer than 7 days may cause localised irritation and should be snipped off or removed as indicated. Avoid prolonged contact of any suture with salt solutions (urinary or biliary), as it may result in calculus formation. Contaminated or infected wounds should be managed with acceptable surgical practice.

An immobilisation of joints by an external support should be considered by the surgeon in the closure of the sites which may require additional support (e.g. expansion, distension, etc.).

Open the outer bag of the foil pack by pulling its two loose ends away from each other - in such a way that the sterile bag containing the suture material falls out undamaged, onto a sterile surface. The sterile bag should then only be touched with sterile gloves and opened by tearing the upper part of the bag.

The product is intended for single use!

When suturing, remember:

- Always double check that you have tied a tight knot. Additional throws should be considered in some surgical circumstances.
- Users should exercise caution when handling surgical needles to avoid inadvertent needle stick injury, which may for example result in transmission of blood-borne pathogens.
- When suturing be careful to avoid damage (suture damage, needle damage, etc.) with surgical instruments (e.g. needle holder, forceps, etc.) or needle. Damage of the suture material could result in negative changes in the tensile strength of the suture material and alter its safety.
- Hold the needle in the area between the third and the half of the distance between the extreme of the union with the thread and the tip. If you hold it near to the tip, you could damage the functional integrity of the needle or break it. If you hold it by the union zone with the thread, the needle body could be broken. When the needle function is altered, it could generate a decrease in its folding resistance and break.
- Discard any unused suture material remaining. Suture needles should be discarded into separated "sharps" container. Dispose of contaminated devices and packaging materials utilizing standard hospital procedures and universal precautions for bio hazardous waste. Do not resterilise!

Do not use if the package is opened or damaged.

This device was designed, tested and manufactured for single patient use only. Reuse, use of the device with opened or damaged packaging, reprocessing and/or resterilisation of this device may lead to fist failure and subsequent injury, illness or death of the patient and/ or create the risk of contamination and patient infection, illness or death of the patient.

Sterilisation

These sutures are sterilised by ethylene oxide or gamma radiation as indicated on the label on the box.

Storage

The product should be protected from direct sunlight and heat, and stored in its original packaging in a clean, dry room at a temperature from 5°C to 30°C. Do not use after

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Müller & Weygandt GmbH · Reichardsweide 40 · 63654 Büdingen/Germany · Tel.: +49 (0) 60 42 - 88 00 88
LUX-SUTURES S.A. · Gruuss-Strooss 22 · 9991 Weiswampach/Luxembourg