Mini LED III

SuperCharged

User's manual
Manuel d'utilisation
Manual de uso
Benutzerhandbuch
Manuale d'uso
Handleiding voor de gebruiker
Manual de utilização
Användningshandbok
Εγχειρίδιο χρήσης

SATELEC
ACTEON


Fig. 1

Fig. 2

| | |
|-------------------|------------|
| ENGLISH | 5 |
| FRANÇAIS | 19 |
| ESPAÑOL | 33 |
| DEUTSCH | 47 |
| ITALIANO | 61 |
| NEDERLANDS | 75 |
| PORTUGUÊS | 91 |
| SVENSKA | 105 |
| ΕΛΛΗΝΙΚΑ | 120 |

ENGLISH SUMMARY

| | |
|-----------------------------------|----|
| I - INTRODUCTION | 6 |
| II - WARNINGS | 6 |
| III - DESCRIPTION | 8 |
| IV - INSTALLATION/STARTUP | 10 |
| V - ROUTINE OPERATION | 11 |
| VI - CONFIGURATION | 11 |
| VII - INTERLOCKS | 11 |
| VIII - MAINTENANCE | 12 |
| IX - TROUBLESHOOTING | 13 |
| X - ELECTROMAGNETIC COMPATIBILITY | 14 |
| XI - DISPOSAL AND RECYCLING | 17 |
| XII - LIABILITY | 17 |
| XIII - ACCESSORIES | 17 |
| XIV - REGULATIONS | 18 |
| XV - SYMBOLS AND ABBREVIATIONS | 18 |

I - INTRODUCTION

You have just taken possession of your Mini LED SuperCharged unit, congratulations.


Created by SATELEC®, the Mini LED SuperCharged is a photopolymerizer allowing dental care to be given using the suitable light guide supplied with the unit.

To take full advantage of the technology of this product, we request that you read carefully the chapter covering all safety recommendations.

The manufacturer's warranty is applicable only if these indications relating to the unit's operation and safety have been correctly applied. All of these safety measures require a sound knowledge of dentistry, photopolymerization and of the specific instructions regarding the operation of the Mini LED SuperCharged given in this operating manual.

Sections with the symbol  points to which we would like to attract your attention.

II - WARNINGS

 **United States Federal Law restricts on its territory the use of this unit exclusively to trained, capable and qualified dental healthcare professionals, or under their control.**

To reduce the risk of accidents, it is imperative to comply with the following precautions:

Users of the unit

Use of the Mini LED SuperCharged is limited exclusively to trained dental healthcare professionals, apt and qualified in the usual context of their activity.

If you have received this unit in error, contact the unit's supplier to arrange for its return.

Interactions/contraindications

Do not use on persons currently suffering from, or that have suffered from in the past, photo-biological reactions (including solar urticaria or erythropoietic protoporphyria) or on persons undergoing treatment using photosensitizing medications (including methoxsalen or chlortetracycline).

Any person, practitioners or patients who have previously suffered from a retina or lens condition or who have undergone eye surgery, in particular cataract surgery, must consult their ophthalmologist before using this unit.

Even in the event of agreement, it is strongly recommended to be prudent as the light intensity may

cause accidents.

It is particularly recommended to continuously wear protective glasses suited to the use of equipment emitting radiation of wavelengths less than 500 nm.

The light radiation produced by this type of equipment can be dangerous and must never be directed towards the eyes, even if the practitioner or the patient is wearing protective glasses suited to the use of equipment emitting radiation of wavelengths less than 500 nm.

The lighting produced by this lamp must be directed only to the part to be treated in the oral cavity.

Interference may occur when used on patients with cardiac pacemakers. This system emits electromagnetic fields, which means there are some potential risks. The malfunctioning of implantable devices such as cardiac pacemakers and ICDs - implantable cardioverter defibrillator (in general, any type of active implant) is possible:

- Ask patients and users if they have an implanted device before using this product. Explain the circumstances to them
- Weigh the risks and benefits and contact your patient's cardiologist or appropriate qualified healthcare professional prior to performing the treatment
- Keep this product away from implanted devices
- Make appropriate emergency provisions and take immediate action if patients become ill
- Symptom including a raised heartbeat, irregular pulse, and dizziness may signal problems with a cardiac pacemaker or ICD.

The unit complies with current electromagnetic compatibility standards, nevertheless, the user shall ensure that any electromagnetic interference does not create an additional risk (presence of radiofrequency emitters, electronic equipment...).

Overexposure to light radiation of the pulp and soft tissues can result in the release of heat and can result in injury to the patient.

To use your equipment in the best possible conditions, it is important to comply with the specifications given in sections II - WARNINGS Equipment users and V - ROUTINE USE.

As far as possible, avoid the accumulation of heat due to a dental dam.

To avoid a sensation of heat, it is recommended to separate cycles by 10 seconds for polymerization times of 10 seconds or less, and by 30 seconds for

polymerization times greater than 10 seconds.

Connection to AC power supply

Ensure that your equipment is connected to the AC power supply by an authorized dentistry installation technician.

Do not tilt the unit at an angle of greater than 5°.

Before connecting the unit, check that the AC voltage is compatible with that indicated on the AC power adapter allowing lamp charging.

A different voltage would result in damage to the unit and could injure the patient and/or the user.

The AC power supply used to power the unit must comply with the applicable standards of your country.

Any fluctuation of the voltage of the ac power supply or electromagnetic field, noncompliant with the applicable limits, may affect the unit's operation.

The device is intended to be connected to a separate power supply (FRIWO FW7660M/12), which is seen as part of the device.

The plug of the power supply is intended to serve as disconnected device, the socket-outlet shall be installed near the equipment and shall be easily accessible.

Unit operation

Do not use the unit if it appears to be damaged or defective.

Do not use the unit if the light guide is damaged (injury hazard...).

Before each use, check that light intensity is compliant (see section 3.1 TECHNICAL DESCRIPTION) using the power verification window on the charging base.

When handling the AC power adapter and/or the battery disconnected from the handpiece, avoid contact with the patients or other persons.

Do not touch the accessible battery and charging base connectors (bump contacts).

To disconnect the ac power supply adapter, grip the AC adapter plug and hold the wall socket.

During dental care, the Mini LED SuperCharged must not be connected to its charging base.

In the event that the unit is not in use or stored, or in the case of a prolonged absence, disconnect the AC adapter from the ac power supply and remove the battery from the body of the unit to protect it from a slow and detrimental discharge.

Do not exert excessive pressure on the unit's LCD screen.

Do not change the battery during use.

Do not short-circuit the battery.

Do not short circuit the charging bump contacts on the charging base.

Do not incinerate the battery (risk of explosion).

Environment

Do not immerse and do not use outside.

Do not place the unit close to a heat source.

The use of solvents, detergents or flammable substances can result in damage or even short-circuits.

Ensure that the power cord connecting the AC adapter to the charging base does not prevent persons from moving freely.

The unit must be stored in its original packaging, in an appropriate place, without danger for persons.

For unit transportation, unscrew the battery and protect the light guide from any unexpected shocks.

Any condensation inside electrical equipment can be dangerous.

If the lamp must be transported from a cold place to a warm place, it must not be used immediately, but only after reaching the ambient temperature.

Never insert or try to insert metallic objects into the equipment to avoid an electric shock or short-circuit hazard.

The unit is not designed for operation in the presence of anesthetic gas or any other flammable gas.

Do not expose the unit to water mist or water splashing.

The unit is not designed for operation close to ionizing radiation.

Maintenance

Before and after each use, it is essential that the unit be disinfected with the products recommended by SATELEC.

Before each use, it is essential that a cleaned and disinfected rigid protection shield be used.

Before each intervention, it is essential that a cleaned, disinfected and sterilized light guide be used.

Before each intervention, check the integrity of the unit and its accessories.

Accessories

Do not use accessories other than those supplied by SATELEC.

The manufacturer refuses to accept any responsibility if damaged parts or accessories are not replaced exclusively by those supplied by the manufacturer. In particular, the use of light guides, ac adapters or batteries other than those supplied by the manufacturer may be dangerous for the patient and the user.

Repair

Do not perform equipment repairs or modifications without the prior authorization of SATELEC.

In the event of an anomaly, immediately disconnect the unit's base and ensure that nobody can use the lamp before verification by the manufacturer or the supplier. This anomaly may be due to noncompliance with safety rules or because of technical damage to the unit.

In the event of an anomaly, contact the supplier of the unit rather than just any repairer, which may return your unit in a hazardous state for both you and your patients.

www.acteongroup.com

Email: satelec@acteongroup.com

III - DESCRIPTION

3.1 PHYSICAL DESCRIPTION

The Mini LED SuperCharged includes the following components:

- Handpiece (fig.1-1).
- Lithium-ion battery (fig.1-8).
- Multi-fiber light guide, 45° curve, Ø 7.5 mm, sterilizable (fig.2-2).
- Charging base with handpiece holder (Fig. 1-9 and 11).
- AC adapter (FRIWO FW7660M/12) (fig.1-12).
- The main plugs (USA, AUS, EEC, UK).
- The rigid protection shield.
- Accompanying documentation.

On option:


- Multi-fiber opalescent light guide, 45° curve, Ø 5.5 mm, sterilizable.
- kit of 5 Flexible protective shields, sterilizable.

3.2 TECHNICAL DESCRIPTION

The Mini LED SuperCharged is equipped with a Light-Emitting Diode (LED) that emits blue visible light of a wavelength between 420 and 480 nm for the photopolymerization of dental materials.

- LCD (Liquid Cristal Display) screen (Fig.2-5)

The LCD screen is backlit and allows the display of the different pieces of information required by the user.

The LCD comprises, from top to bottom:

- Display of time counter for polymerization light emission time, expressed in seconds.
- A pictogram indicating the Fast, Pulse, and Stepped menus.

| Pictogram | Modes |
|-----------|-----------------|
| | Fast-Cure mode |
| | Pulse-Cure mode |
| | Step-Cure mode |

- Display of the battery charge level status (battery being most charged when all of the pictogram's round symbols are displayed and the least charged when all of the round symbols have disappeared from the LCD screen).

- Control buttons (fig. 2 - 3/6/7)

The handpiece comprises three buttons:

- An ON/OFF button (fig.2-3) to start or stop the polymerization cycle.
- Left button (fig.2-6) for selecting modes.
- Right button (fig. 2-7) for selecting the exposure time.

- Types of modes

You can choose three types of modes by pushing the left button (fig. 2-6):

- Exposure times

Once the mode has been selected, you can choose the exposure time by pushing the right button (fig.2-7):

- Fast-Cure mode (full power):
 - 3 seconds (show 03) ;
 - 4 seconds (show 04) ;
 - 5 seconds (show 05) ;
 - 10 seconds (show 10) ;
- Pulse-Cure mode:
 - 5 shots of 1 second (full power with emission of 5 successive one second flashes with a rest of period of 250 ms between flashes) - (show 05);
 - 10 shots of 1 second (full power with emission of 10 successive one second flashes with a rest of period of 250 ms between flashes) - (show 10) ;
- Step-Cure mode:
 - 6 seconds progressively and 3 seconds at full power (show 09) ;
 - 10 seconds progressively and 5 seconds at full power (show 15).

- Types of warnings

An audible alarm (beep) is triggered at the start of polymerization.

- Charge time

Battery charge time is about **3 hours**.

- Battery charge level indicator

When only one element of the pictogram remains (the one farthest to the right on the LCD screen), two beeps warn you that there are about 500 seconds of polymerization remaining.

The operator is then invited to return the Mini LED SuperCharged to the charging base.

When the battery is completely discharged, the LCD screen displays the message "Lb" (Low battery) and no round symbols are shown, the handpiece then beeps 4 times and the LED indicator light becomes red.

- Indicator light

The handpiece comprises 1 two-color indicator light (fig.2-4) located on the ON/OFF button.

When the unit detects a fault (see section IX Troubleshooting), the indicator light becomes continuously red. In all other cases, it remains continuously green (except when the handpiece is placed in standby mode, in which case the indicator light goes out).

- Charging base

The charging base, including the translucent handpiece holder (fig.1-9), which also serves as a verification indicator, is used to recharge the battery.

Place the charging base on a hard and stable surface, with an inclination of no more than 5°.

The translucent handpiece holder base (fig.1-9) flashes three times (red, green, blue) and beeps when connected to ac power.

Place the lamp on the holder ensuring that there is good contact between the bump contacts and the accessible connectors under the battery.

Two beeps confirm that the lamp is correctly positioned on its holder.

The translucent holder base lights up, with a blue color, and starts flashing, indicating that the battery is charging.

The battery is completely charged when the translucent holder base stops flashing and the blue light remains lit continuously.

A power level verification window is integrated in the charging base.

This power level verification window allows you to precisely check if your lamp, equipped with its light guide, is sufficiently powerful.

The input window of the verification system is located on the front of the holder (fig.2-10).

- Power level verification window

The procedure for using the power level verification window is as follows:

Program the lamp for 3-second Fast mode.

Insert the light guide into the handpiece after having checked that it is intact and clean.

In the same way, first check that the power level verification window is intact and clean.

Place the extremity of the light guide flat against the power level verification window (fig.2-10) and activate your Mini LED SuperCharged.

According to the power level detected, the translucent handpiece holder base (fig.2-9) will become green in the case of a compliant power level or red if the power is insufficient.

In the event of the power being insufficient, please refer to chapter IX TROUBLESHOOTING.

- technical specifications

| | |
|--|---|
| Model name: | Mini LED SuperCharged |
| Medical category: | Ila according to the european directive 93/42/CEE |
| Dimensions without light guide: | |
| Weight: | 160 g |
| Dimensions: | Ø24 x 201 mm |
| Operation: | Continuous operation |
| Protection: | |
| Category: | Type B |
| Protection: | 5 AT FU1 fuse (non accessible) 125 V |
| Protection index: | IPX0 |
| AC power adapter: | |
| Input voltage: | 100 V AC to 240 V AC |
| Frequency: | 50 Hz to 60 Hz |
| Output voltage: | 12 V DC |
| Output current: | 0.8 A |
| Category: | II |
| Protection index: | IP 41 |
| Reference | FRIWO FW7660M/12 |

Charging base:

Input voltage: 12 V DC
 Protection: 3 AT F1 fuse (non accessible) 125 V
 Category: Continuous operation
 Protection index: IPX0

Battery:

Type: Lithium-Ion
 Dimensions: 88 x Ø24 mm
 Capacity: 2300 mAh

Optical specifications:

- LED for polymerization:
 Wavelength range: 420 - 480 nm
 Central wavelength: 455 - 465 nm
 Intensity: 2000 mW/cm² ± 10% for an active fiber diameter of 7.5 mm
 Maximum exposure time: 15 seconds (show 15)

Temperatures:

Operation: +10°C to +40°C
 Storage: -20°C to +70°C

Humidity:

Operation: 30% to 75%
 Storage: 10% to 100%
 condensation included

Atmospheric pressure 700hPa to 1060hPa

IV - INSTALLATION/STARTUP

4.1 UNPACKING THE UNIT

On reception of the unit, look for any damage that may have occurred during transportation. If necessary, contact your supplier.

4.2 RECOMMENDATIONS

Check that the environmental conditions have been complied with (ambient temperature between 10°C and 40°C and humidity between 30% and 75%).

4.3 INSTALLATION

 **Ensure that the unit is not installed adjacent to or on top of any other equipment. Do not place the power cord in a cable feedthrough or cable cover.**

Remove the protective caps from the handpiece, screw the battery onto the handpiece and insert the sterilized light guide into the handpiece. Ensure that the light guide is properly inserted, confirmed by a click. Ensure that all the segments of the LCD screen are present once the battery has been properly screwed down. Place the base on a hard and stable surface that is not inclined by more than 5°. Connect the AC adapter after ensuring that the voltage indicated corresponds to the electrical installation available in the dental practice. Two audible beeps confirm that the lamp has been correctly placed on its holder. The translucent holder base lights up in blue and starts flashing, indicating that the battery is charging. As soon as the battery is charged, the translucent holder base stops flashing and the blue light remains illuminated continuously.

4.4 OPERATING FOR THE FIRST TIME

Before using for the first time:
 - The delivered batteries are loaded at 60% before shipping. It is therefore necessary to proceed to a complete loading.
 - it is essential that the sterilizable accessories (light guide and protective shield) be sterilized and that the unit base be disinfected (see chapter concerning maintenance).

4.5 MODES AVAILABLE

The Mini LED SuperCharged has three modes you can choose by pushing the left button (Fig.2-6) :

- Fast-Cure,
- Pulse-Cure,
- Step-Cure.

4.6 EXPOSURE TIMES AVAILABLE

In these modes, you can choose different exposure times. By pushing the right button (fig.2-7), you can select the following combinations:

- Fast-Cure mode:
 - 3 seconds (show 03) ;
 - 4 seconds (show 04) ;
 - 5 seconds (show 05) ;
 - 10 seconds (show 10) ;
- Pulse-Cure mode:
 - 5 shots of 1 second (full power with emission of 5

- successive one second flashes with a rest of period of 250 ms between each flash) - (show 05);
- 10 shots of 1 second (full power with emission of 10 successive one second flashes with a rest of period of 250 ms between each flash) - (show 10) ;
- Step-Cure mode:
 - 6 seconds progressively and 3 seconds at full power (show 09) ;
 - 10 seconds progressively and 5 seconds at full power (show 15).

V - ROUTINE OPERATION

The Mini LED SuperCharged is normally placed on its holder. The first time you plug it into the battery it will activate in the Menu Fast-Cure and the exposure time of 3 seconds.

Your Mini LED SuperCharged is ready to operate once the menu and the exposure time has been selected. Place the light guide as close as possible to the surface of the material to be polymerized, without touching it as this may adversely affect the quality of the polymerization. The use of a transparent strip is recommended. The polymerization starts by pushing on the ON/OFF button. The setting will be confirmed by an audible signal (beep). A countdown on the LCD screen shows you how much time remains. When the polymerization cycle is completed, the time that was applied is displayed. You can press the ON/OFF button at any time to interrupt the current polymerization cycle. After 5 minutes without use, the unit goes into standby mode (low consumption). The green indicator light and the LCD screen backlighting switch off. The unit can be taken out of standby mode by simply pressing one of the three buttons (which will not activate that function when coming out of standby mode).

VI - CONFIGURATION

6.1 STARTUP

On startup, the unit performs an automatic test sequence (auto-check).

The device automatically recalls the settings used for the previous dental treatment.

6.2 MODES

The Fast, Pulse and Stepped modes are selected using the left button (fig.2-6) on the handpiece.

6.3 EXPOSURE TIMES

The exposure times are selected by pressing the right button (fig.2-7).

VII - INTERLOCKS

The Mini LED SuperCharged is equipped with a system for the detection of possible unit operation anomalies.

7.1 EXCESSIVE TEMPERATURE

During intensive use, an excessive temperature may be detected and the screen then displays the letters "OH" (Over Heat), the audible warning (beep) sounds 4 times and the indicator light (two-color LED) becomes red.


It is then recommended that the user leaves the handpiece to cool for a few minutes, until the indicator light becomes green again and the screen once again displays the selected time cycle.

7.2 BATTERY CHARGE LEVEL

When the LCD screen displays the last remaining round battery charge indicator, the audible alarm (beep) sounds twice.

You then have about 500 seconds remaining before the unit becomes completely discharged. It is then recommended, if possible, to charge the battery for further use.

When the battery becomes completely discharged, the letters "Lb" (Low Battery) are displayed on the screen, the audible alarm (beep) sounds 4 times and the indicator light (two-color LED) becomes red.

The user must return the handpiece to its charging base to recharge the unit's battery (see chapter 4.3 INSTALLATION).

VIII - MAINTENANCE

 Before conducting any maintenance on the Mini LED SuperCharged, check that:

- The Mini LED SuperCharged is not on its charging base.
- The battery has been unscrewed from the handpiece.
- The charging base has been disconnected from the AC supply.

Before cleaning the handpiece, insert the protective cap in the place of the light guide, supplied to ensure that liquid does not enter the handpiece.

Avoid using cleaning and disinfection products containing flammable agents (or other corrosive agents such as acetone, chlorine or bleach). Otherwise, ensure that the product completely evaporates and that there are no combustibles on the unit and its accessories before operation.

Do not use abrasive products to clean the unit.

Do not immerse the unit.

Never use ultrasound cleaning, whether for the Mini LED SuperCharged or its accessories.

None of the accessories are delivered in a sterilized state.

Only the light guide and the flexible protective shield are sterilizable.

Before sterilization, check the cleanliness of your autoclave and the quality of the water used.

After each sterilization cycle, immediately remove the items from the autoclave to reduce the risk of corrosion of metallic parts.

It is necessary to leave the sterilized items to cool down to the ambient temperature and dry before re-using them.

It is recommended to sterilize the items identified as sterilizable by Satelec individually in the sterilization bags foreseen for this purpose.

To maintain the sterile and aseptic state of accessories, make sure that they are kept in hermetically-sealed bags or containers suited to use

in dentistry.

The maintenance and/or sterilization instructions that follow must be applied before each use of the unit.

8.1 - Pre-disinfection/cleaning

Clean and disinfect the body, light guide, flexible shield and rigid protection shield, and the charging base of the Mini LED SuperCharged using ready-to-use cleaning/disinfecting wipes based on alcohol, amphoteric disinfectant and biguanide (refer to the manufacturer's instructions) for at least two minutes.

Leave the product to act for at least 15 minutes.

Use wipes with CE marking, or compliant with any standard that may be required by national regulations.

8.2 - Drying

Dry using a clean, single-use, non-woven cloth to remove any liquid traces.

8.3 - Packing

Pack in single-use sterilization bags or sleeves compliant with the requirements defined in the ISO 11607-1 standard, or compliant with any standard that may be required by national regulations.

8.4 - Sterilization

The optic fiber and the flexible protective shield for the Mini LED SuperCharged must be sterilized individually in an autoclave according to the following parameters:

- Autoclave: Type B compliant with the EN 13060 standard.
- Sterilization temperature: 134°C.
- Sterilization steady state: 18 minutes.
- Pressure: 2 bar minimum.

CAUTION: Not all autoclaves can reach 134°C. Not all autoclaves perform a pre-depressurization. For further information on the applicable sterilization instructions, consult the autoclave manufacturer.

8.5 - Storage

After this, store the sterilized items in a dry, dust-free place.

Before re-use, in the event of nonconforming packaging integrity, re-package and re-sterilize according to the defined protocol.

IX - TROUBLESHOOTING

In the event of a problem, before contacting the after-sales service of SATELEC or the supplier:

- Ensure that the base is correctly connected to the AC supply to ensure that the battery is charged normally. If the Mini LED SuperCharged s holder is not lit despite being correctly connected, contact the after-sales service of SATELEC or the supplier.
- Check that at least one of the four round battery charge level indicators is visible on the LCD screen before pressing the ON/OFF button.
- In the event of a faulty battery (see VII - INTERLOCKS) an automatic interlock system will prevent the unit from operating. If this is the case, recharge the battery by placing the lamp on its holder or by using a second fully-charged battery.
- Intensive use of the lamp can result in high temperatures inside the unit. If this happens, an automatic interlock system prevents the unit from operating (see VII - INTERLOCKS). Leave the unit at rest for a few minutes to allow it to cool down.
- After each use, check that there is no composite residue adhering to the light guide. If this is the case, immediately remove the residues and check that the surface of the light guide has not been damaged. If damage is visible, replace the light guide as the lamp's power could be significantly reduced.
- Under normal conditions of use, the lamp's power does not vary if the battery is properly charged. Consequently, it is not necessary to check power as with ordinary polymerization lamps. Nevertheless, in case of doubt, check its power using your light guide.
- The battery located at the bottom of the handpiece can be replaced by unscrewing it.
- After removing the light guide, check that the LED is clean and undamaged. If necessary, clean it using a medical-quality dry air jet (free from compressor oil).
- During verification of the power level of the handpiece, if the translucent handpiece holder base is red, check if the LED, light guide and verification window are clean.
- In the event of dust being present, clean using a dry air jet.
- If the problem persists or if the light guide or the verification window are damaged, the lamp must be returned to the after-sales service.

 In the event of an anomaly, contact the supplier of the unit rather than just any repairer, which may return your unit in a hazardous state for both you and your patients.

The technical service of your supplier is at your disposal for any technical problems encountered on your unit.

| Anomaly observed | Possible causes | Solutions |
|--|--|--|
| No operation (LCD screen off) | Battery completely discharged | Recharge battery |
| | Battery defective | Return to SATELEC after-sales department |
| | MiniLED defective | Return to SATELEC after-sales department |
| No operation (LCD screen on) | MiniLED defective | Return to SATELEC after-sales department |
| Light power defect or no light power | LED defective and/or dirty | Check cleanliness of verification window and/or Return to SATELEC after-sales department |
| | MiniLED defective | Return to SATELEC after-sales department |
| | Defective buttons | Return to SATELEC after-sales department |
| | Light guide defective and/or dirty | Clean light guide and/or Return to SATELEC after-sales department |
| | Reflector defective or dirty | Clean reflector (dry air jet) and/or Return to SATELEC after-sales department |
| Charger does not operate | AC power wall outlet defective | Contact your electrician |
| | AC power adapter defective | Return to SATELEC after-sales department |
| | Light and/or audible alarms defective | Return to SATELEC after-sales department |
| | Contact failure on Jack connector socket | Return to SATELEC after-sales department |
| | Fuse defective | Return to SATELEC after-sales department |
| Dysfunction of power level verification | Indicator light defective | Return to SATELEC after-sales department |
| | Light guide defective and/or dirty | Clean light guide and/or Return to SATELEC after-sales department |
| | LED defective and/or dirty | Clean LED (dry air jet) and/or Return to SATELEC after-sales department |
| | Reflector defective and/or dirty | Clean reflector (dry air jet) and/or Return to SATELEC after-sales department |
| | Window defective and/or dirty | Clean window and/or Return to SATELEC after-sales department |

X - ELECTROMAGNETIC COMPATIBILITY

Warning: The charger power cord must be kept apart from those of any nearby devices.

Mini LED SuperCharged requires special precautions to be taken with regard to electromagnetic compatibility. It must be installed and prepared for use as described in chapter 3.

Certain types of mobile telecommunication devices such as mobile telephones are likely to interfere with the Mini LED SuperCharged. The recommended separation distances in this paragraph must therefore be complied with.

The Mini LED SuperCharged must not be used near or on top of another device.
 If this cannot be avoided, its operation under the conditions of use must be checked beforehand.
 The use of accessories other than those specified or sold by SATELEC as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the Mini LED SuperCharged.

10.1 - Electromagnetic emissions

The Mini LED SuperCharged is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the Mini LED SuperCharged is used in such an environment.

| Emission test | Compliance | Electromagnetic environment - guidance |
|-------------------------|--------------------|---|
| RF emission - CISPR 11. | Group 1 Class B | Mini LED SuperCharged uses RF energy for internal operation. Therefore, its radiofrequency emissions are very low and are not likely to cause any interference in nearby equipment. Mini LED SuperCharged is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes. |

10.2 - Magnetic and electromagnetic immunity

The Mini LED SuperCharged is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the Mini LED SuperCharged is used in such an electromagnetic environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---|--|--|---|
| Magnetic field at 50Hz. IEC 61000-4-8 | 3A/m | 3A/m | The intensity of the magnetic field should be equivalent to that of a typical commercial or hospital environment (hospital, clinic). |
| Electrostatic discharge (ESD) IEC 61000-4-2. | ± 6KV contact ± 8KV air | ± 6KV contact ± 8KV air | Floors must be wood, concrete, cement or tiled. If floors are covered with synthetic material (carpet, etc.), the relative humidity must be at least 30%. |
| Electrical fast transients IEC 61000-4-4. | ± 2KV for power supply lines | ± 2KV for power supply lines | Power quality should be that of a typical commercial or hospital environment. |
| Surges IEC 61000-4-5. | ± 1KV differential mode ± 2KV common mode | ± 1KV differential mode ± 2KV common mode | Power quality should be that of a typical commercial or hospital environment. |
| Voltage dips, short interruptions and voltage variations IEC 61000-4-11. | <5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles | <5% UT (>95% dip in UT) for 0.5 cycles. 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 250 cycles | Power quality should be that of a typical commercial or hospital environment. If the use of the Mini LED SuperCharged requires continued operation during a power cut, it is recommended that the product be powered from a separate power supply (UPS, etc.). |

10.3 - Electromagnetic immunity / mobile radiofrequency equipment

The Mini LED SuperCharged is intended for use in the electromagnetic environment specified in the table below. The user and/or installer must ensure that the device is used in such an electromagnetic environment.

| Immunity test | IEC 60601 test level | Compliance level | Electromagnetic environment - guidance |
|---|----------------------------|------------------|---|
| Portable and mobile radiofrequency communications devices must not be used near the Mini LED SuperCharged (including its cables) at a distance less than that recommended and calculated according to the frequency and power of the emitter. | | | |
| Conducted disturbance, radiofrequency fields. IEC61000-4-6 | 3 V/m 150 KHz to 80 MHz | 3 V/m | Recommended separation distance: $d = 1.2 \sqrt{P}$ |
| Radiated radiofrequency electromagnetic field. IEC61000-4-3 | 3 V/m 80 MHz to 2.5 GHz | 3 V/m | $d = 1.2 \sqrt{P}$ 80MHz to 800MHz. $d = 2.3 \sqrt{P}$ 800MHz to 2.5GHz. Where P is the maximum power rating of the emitter in watts (W) according to the manufacturer's specifications and d is the recommended minimum separation distance in meters (m). |
| The electromagnetic field strengths of fixed radiofrequency emitters, as determined by an electromagnetic environment measurement (a), must be less than the compliance level in each frequency range (b). Interference may occur near equipment marked with the symbol below:  | | | |

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

(a): The electromagnetic field strengths of fixed radiofrequency emitters, such as base stations for mobile telephones (cellular/cordless), mobile radios, amateur radios, AM/FM radio broadcasts and TV broadcasts cannot be determined exactly by theory. To assess the electromagnetic environment due to fixed radiofrequency emitters, an electromagnetic environment measurement must be made. If the measured radiofrequency field strength in the immediate environment where the product is used exceeds the compliance level specified above, the performance of the product must be tested to verify whether it conforms to the specifications. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the product.

(b): In the 150 kHz to 80 MHz frequency range, the electromagnetic field strengths must be less than 3 V/m.

10.4 - Recommended separation distances

The Mini LED SuperCharged is intended for use in an electromagnetic environment in which radiated radiofrequency disturbances are controlled.

The Mini LED SuperCharged user and/or installer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile radiofrequency communications equipment (emitters) and the Mini LED SuperCharged, according to the maximum output power of the equipment, as recommended in the table below.

| Rated max. power of the emitter (W) | Separation distance in meters (m) according to emitter frequency | | |
|-------------------------------------|--|---|--|
| | 150 KHz to 80 MHz $d = 1.2 \sqrt{P}$ | 80 MHz to 800 MHz $d = 1.2 \sqrt{P}$ | 800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$ |
| 0.01 | 0.12 m | 0.12 m | 0.23 m |
| 0.1 | 0.38 m | 0.38 m | 0.73 m |
| 1 | 1.2 m | 1.2 m | 2.3 m |
| 10 | 3.8 m | 3.8 m | 7.3 m |
| 100 | 12 m | 12 m | 23 m |

For emitters rated at max. power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the emitter, where P is the max. power rating of the emitter in watts (W) according to the manufacturer.

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These specifications may not be applicable in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and persons.

10.5 - Cable lengths

| Cables and accessories | Maximum length | Complies with: |
|------------------------|----------------|--|
| charger power cord | <3m | RF emission, CISPR 1 - Class B Immunity to magnetic fields: IEC61000-4-8. Immunity to electrostatic discharge: IEC61000-4-2 Immunity to electrical fast transients/bursts: IEC61000-4-4 Immunity to surges: IEC61000-4-5 Immunity to voltage dips, short interruptions and voltage variations: IEC61000-4-11 Immunity to conducted disturbances induced by radiofrequency fields: IEC61000-4-6 Immunity to radiated radiofrequency electromagnetic fields: IEC61000-4-3 |

XI - DISPOSAL AND RECYCLING

As electrical and electronic equipment, the device must be disposed of according to a specialized procedure for collection, pick-up and recycling or destruction (in particular on the European market, with reference to Directive 2002/96/EC of 27/01/2003).

When your device reaches the end of its life, we consequently recommend that you contact your dental equipment dealer (or, failing this, the nearest ACTEON GROUP office, the list of which is given in chapter 16), for information on how to proceed.

XII - LIABILITY

The manufacturer is not liable if:

- the manufacturer's installation recommendations have not been followed (supply voltage, electromagnetic environment, etc.);
- repairs have been performed by persons not authorized by the manufacturer;
- the device has been used in an electrical installation which does not comply with current standards;
- the device has been used in a way which is not stipulated in this Manual;
- accessories other than those supplied by SATELEC have been used;
- the instructions in this document have not been followed.

The manufacturer reserves the right to modify the unit and/or the Operating Manual without notice.

XIII - ACCESSORIES

The following accessories are available for the Mini LED SuperCharged:

- Opalescent light guide Ø 7.5 mm : Ref. F02648
- Opalescent light guide Ø 5.5 mm(option) : Ref. F02652
- Flexible protective shield x 5 (option): Ref. F61507
- Rigid protection shield: Ref. F05407
- AC power supply: Ref. F05216
- Handpiece: Ref. F05215
- Battery: Ref. F05211


XIV - REGULATIONS

This medical device is classified as class IIa according to European Directive 93/42/EEC.

This equipment is manufactured in compliance with the current IEC 60601-1 standard.

This equipment has been designed and manufactured according to an ISO 13485-certified quality assurance system.

XV - SYMBOLS AND ABBREVIATIONS

| SYMBOL | DEFINITION |
|---|--|
|  | Alternating current |
|  | Direct current |
|  | Do not discard as household waste |
|  | Follow operating instructions |
|  | "ON" / "OFF" (pushbutton) |
|  | Type B |
|  | Protective glasses recommended |
|  | CE marking |
|  | No access for persons with active implant as pacemaker |

Note:

Technical personnel of the Satelec authorized dealer network can obtain from **ACTEON Group** on request all the information they need for repair of the parts of the curing light that Satelec has identified as repairable.

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